



SECURESAFE™ SAFETY PEN NEEDLE

Dynarex SecureSafe Safety Pen Needles are designed to deliver injectable medications safely and efficiently. A retractable shield helps reduce the risk of needlestick injuries. Compatible with most insulin pens, they offer added convenience in clinical and home settings.

- 9655: Dual-safe design with automatic locking safety mechanisms at patient end and pen-hub end helps reduce risk of needlestick injuries
- 9656: Retractable safety shield helps reduce the risk of needlestick injuries
- Tri-bevel cut tip supports gentle and effective insertion
- Film-coated surface helps promote a smoother injection
- Compatible with most insulin pens for optimal versatility
- Sterile, non-toxic, and non-pyrogenic
- Not made with natural rubber latex
- 100 per box, 12 boxes per case

NEW	9655	31G × ³ / ₁₆ " (0.25 MM × 5 MM)
	9656	30G × ⁵ / ₁₆ " (0.31 MM × 8 MM)



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |