

GLOVE COMPARISON CHART



SafeTouch™

LATEX GLOVES	Item	Color	Thickness	Textured	Tensile Strength	Meets Standards	Length	Width	Sterile	Chemo Approved	Low Protein	Non-Chlorinated	Premium Chlorinated	Polymer Lined	NFPA Certified
Examination Gloves	2336-2339		3 mil min.	Fingertip	18 MPa, Min.	ASTM D3578	220 mm minimum for S; 230 mm minimum for M, L, XL	80±10 for S; 95±10 for M; 111±10 for L; 115±10 for XL						●	
Sterile Surgical Gloves	2360-2390		4 mil min.	Palm	26-28 MPa, Min.	ASTM D3577	275 mm minimum for all sizes	76±6 for Size 6; 82±6 for Size 6.5; 88±6 for Size 7; 92±6 for Size 7.5; 101±6 for Size 8; 106±6 for Size 8.5; 112±6 for Size 9	●						
Sterile Examination Gloves	2451-2453		3 mil min.	Palm	19-21 MPa, Min.	ASTM D3578	245 mm minimum for all sizes	85±5 for S; 96±5 for M; 106±5 for L	●						
Black Arrow™ Examination Gloves	2321-2324		6 mil min. at palm; 7 mil min. at finger	Palm	≥21 MPa, Min.	ASTM D3578	230 mm minimum for all sizes	94 for S; 105 for M; 116 for L; 125 for XL						●	
Accutouch® Examination Gloves	6621-6625		4 mil min. at palm; 5 mil min. at finger	Palm	14 MPa, Min.	ASTM D3578	220 mm minimum for XS, S; 230 mm minimum for M, L, XL	70±10 for XS; 80±10 for S; 95±10 for M; 111±10 for L; 115 ±10 for XL				●		●	
Aloetex Examination Gloves	6716-6719		6 mil min. at palm; 7 mil min. at finger	Palm	≥21 MPa, Min.	ASTM D3578	230 mm minimum for all sizes	77±5 for XS; 83±5 for S; 94±5 for M; 105±5 for L						●	

- All our gloves are FDA 21 CFR Part 820 compliant, and registered with the FDA in their Medical Device Listing database.
- Latex gloves are not recommended for use in the food industry, as some people may be allergic to latex material. Therefore, vinyl gloves are used for most food industry applications.
- Dynarex reserves the right to make changes to these statements

GLOVE COMPARISON CHART

SafeTouch™

NITRILE GLOVES	Item	Color	Thickness	Textured	Tensile Strength	Meets Standards	Length	Width	Sterile	Chemo Approved	Low Protein	Non-Chlorinated	Premium Chlorinated	Polymer Lined	NFPA Certified
Blue Nitrile Examination Gloves	2511-2514		3-4 mil at palm; 4-6 mil at finger	Fingertip	14 MPa, Min.	ASTM D6319; ASTM D5151; ASTM D6978-05	240 mm minimum for all sizes	85±10 for S; 98±10 for M; 107±10 for L; 115±10 for XL						●	
Black Nitrile Examination Gloves	2521-2524		3-4 mil at palm; 4-6 mil at finger	Fingertip	14 MPa, Min.	ASTM D6319; ASTM D5151	240 mm minimum for all sizes	85±10 for S; 95±10 for M; 106±10 for L; 115±10 for XL; 125±10 for XXL						●	
DynaPlus Examination Gloves	2515-2519		2 mil min.	Fingertip	≥18 MPa, Min.	ASTM D6319; ASTM D6124	220 mm minimum for XS, S; 230 mm minimum for M, L, XL	70±10 for XS; 85±10 for S; 95±10 for M; 105±10 for L; 115±10 for XL							
Sterile Examination Gloves	6522-6524		2 mil min.	Palm	20-22 MPa, Min.	ASTM D6319; ASTM D5151	240 mm minimum for all sizes	85±10 for S; 97±10 for M; 105±10 for L	●						
Sterile Examination Gloves	6528-6530		2 mil min.	Palm	20-22 MPa, Min.	ASTM D6319; ASTM D5151	240 mm minimum for all sizes	85±10 for S; 97±10 for M; 105±10 for L	●						
Sterile Surgical Gloves	6533-6539		5-8 mil	No	14 MPa, Min.	ASTM D6319; ASTM D5151	280 mm minimum for all sizes	77±5 for Size 6; 83±5 for 6.5; 89±5 for 7; 95±5 for 7.5; 102±6 for 8; 108±6 for 8.5; 114±6 for 9	●					●	
AloeSkin Examination Gloves	6721-6724		2 mil min.	No	14 MPa, Min.	ASTM D6319	220 mm minimum for XS, S; 230 mm minimum for M, L, XL	70±10 for XS; 80±10 for S; 95±10 for M; 110±10 for L; 120±10 for XL						●	

- All our gloves are FDA 21 CFR Part 820 compliant, and registered with the FDA in their Medical Device Listing database.
- Latex gloves are not recommended for use in the food industry, as some people may be allergic to latex material. Therefore, vinyl gloves are used for most food industry applications.
- Dynarex reserves the right to make changes to these statements

GLOVE COMPARISON CHART

SafeTouch™

NITRILE GLOVES	Item	Color	Thickness	Textured	Tensile Strength	Meets Standards	Length	Width	Sterile	Chemo Approved	Low Protein	Non-Chlorinated	Premium Chlorinated	Polymer Lined	NFPA Certified
Gloves in a bag	2618 2619 2624 2625		2 mil min.	Fingertip	≥14 MPa, Min.	ASTM D6319	230 mm minimum	110±10 for L 120±10 for XL							
Gloves in a box	2626UB-2		2 mil min.	Fingertip	≥14 MPa, Min.	ASTM D6319	230 mm minimum	120±10 for X L							
VINYL GLOVES															
Examination Gloves	2611-2614		3 mil min.	No	≥14 MPa, Min.	ASTM D5250	230 mm minimum for all sizes	85±10 for S; 95±10 for M; 105±10 for L; 115±10 for XL						●	
Multi Care Synthetic Examination Gloves	6812 6814 6815		3 mil min.	No	≥14 MPa, Min.	ASTM D5250	230 mm minimum for all sizes	85±10 for S; 95±10 for M; 105±10 for L; 115±10 for XL						●	
Next Generation Stretch Examination Gloves	6824		3 mil min.	No	≥14 MPa, Min.	ASTM D5250	230 mm minimum for all sizes	85±10 for S; 95±10 for M; 105±10 for L; 115±10 for XL						●	
Gloves in a bag	2616-2617		3 mil min. at palm; 2 mil min. at finger	No	≥11 Mpa, Min.	ASTM D5250	230 mm minimum	105±10 for L							

- All our gloves are FDA 21 CFR Part 820 compliant, and registered with the FDA in their Medical Device Listing database.
- Latex gloves are not recommended for use in the food industry, as some people may be allergic to latex material. Therefore, vinyl gloves are used for most food industry applications.
- Dynarex reserves the right to make changes to these statements