



# WOUND, SKIN, AND HYGIENE SOLUTIONS





## ABOUT

---

Dynarex features a comprehensive range of wound care, skin care, and hygiene supplies developed by expert product specialists who understand the challenges providers face while administering proper care. Since our founding in 1967, Dynarex has been one of the industry's leaders in disposable and durable medical equipment. As your partner in the marketplace, we continue to diversify our portfolio across the entire continuum of care because the success of your business matters to us. When it comes to doing business with Dynarex, we choose to work like we're family, because We Care Like Family™.

## LOCATIONS

---

### Corporate Headquarters

155 Chestnut Ridge Road  
Suite 2202  
Montvale, NJ 07645

### International China Office

Tianxiang International Mansion  
No. 318, Suya Road  
Suzhou Industrial Park  
Suzhou, Jiangsu, China

### Dynarex Technologies Pvt Ltd

8th Floor, D-199, Phase 8b  
Industrial Area, Sector 74  
Sahibzada Ajit Singh Nagar  
Punjab 160055 India

### Middletown Warehouse

11 Dynarex Drive  
Middletown, NY 10941

### Arizona Warehouse

7811 N. Glen Harbor Boulevard  
Glendale, AZ 85307

### Arkansas Warehouse

2807 Highway 67 S  
Pocahontas, AR 72455

### Florida Warehouse

14500 Hyatt Road  
Jacksonville, FL 32218

### Ohio Warehouse

702 E. Montford Avenue  
Ada, OH 45810



## ADDITIONAL INFO

Disclaimer: Dynarex will make every effort to inform our customers of any changes made to a specific product, catalog copy or marketing collateral materials. Dynarex cannot be held responsible for these changes. Dynarex reserves the right to make changes to the technical specs, dimensions, look, or any other product change at any time without prior notice except as required by the US FDA Code of Federal Regulations for changes made to medical devices, OTC drugs and other foreign governmental regulations.

Go to [www.dynarex.com](http://www.dynarex.com) for the latest product information. Customers should consult with their Customer Care Representative before placing a purchase order. Please note: All products in the Dynarex catalog may not be licensed in all countries. Customers are hereby directed to consult with their Customer Care Representative before placing a purchase order. Product sizes may vary slightly based on production tolerances using the metric system and are stated for marketing purposes only. Go to [www.dynarex.com](http://www.dynarex.com) to find which products are labeled with the California Proposition 65 warning to indicate it contains chemicals that could cause cancer, birth defects or other reproductive harm. Learn more about prop65 at [www.p65warnings.ca.gov](http://www.p65warnings.ca.gov).

Copyright © 2026 Dynarex Corporation. All rights reserved.

R251010

## ADVANCED WOUND CARE

### WOUND DRESSINGS

CuraFoam Foam Dressings .....	08
DynaFoam Waterproof Bordered Foam Dressing.....	09
DynaLevin Waterproof Adhesive Bordered Foam Dressings.....	09
DynaGuard Waterproof Composite Dressings .....	10
Island Dressings.....	11
FoamFlex Non-Adhesive Waterproof Foam Dressings.....	11
DynaSorb Super Absorbent Dressings.....	12
SiliGentle Silicone Foam Dressings.....	13
SiliGentle AG Silicone Foam Dressings .....	15
DynaFoam AG Bordered Foam Dressing.....	16
CuraFoam AG Foam Dressing.....	17
Oil Emulsion Dressings.....	17
Petroleum Gauze Dressings.....	18
XeroForm Petrolatum Gauze Dressings.....	18
DynaGinate Calcium Alginate Dressings .....	19
DynaGinate AG Calcium Alginate Dressings.....	19
DynaDerm Hydrocolloid Dressings.....	21
Unna Boot Bandages.....	22
View Guard Transparent Dressings.....	23
Packing Strips.....	24
Dressing Retention Sheet.....	25
View Guard Transparent Dressing Rolls.....	25
DynaCol Bovine Collagen Powder.....	26
DynaCol Collagen Wound Dressing.....	26
DynaCol AG Silver Collagen Wound Dressing.....	27

DynaGel Hydrogel Impregnated Gauze Dressing .....	27
Transparent Film Window Dressing.....	28

### MEASURING TOOLS

DynaRule Disposable Paper Rulers .....	28
DynaRule Bullseye Measuring Guide.....	29
Retractable Tape Measure .....	29

### GELS + CLEANSERS

DynaGel Silver Wound Gel Rx.....	30
L-Mesitran Soft.....	31
DynaWound Wound Cleanser Spray .....	31

### NEGATIVE PRESSURE WOUND THERAPY

Negative Pressure Wound Therapy Pump.....	32
Carry Case for NPWT Pump .....	33
Wound Vac Disposable Canisters.....	33
Connectors for Negative Pressure Wound Therapy Device.....	34
Foam Dressing.....	34
Transparent Wound Drape .....	35
Wound Vac Dressing Kit .....	35

### CLINICAL SKIN CARE

Calmara - External Analgesic.....	38
Antifungal Creams.....	39
Vitamins A&D Ointments.....	39
DynaShield Skin Protectant Barrier Cream.....	40

Hydrocortisone Cream USP.....	40
Triple Antibiotic Ointment.....	41
Bacitracin Ointment USP.....	41
AquaSoothe Hydrating Ointment .....	42
Calasoothe Skin Protectant Moisture Barrier Cream .....	42
Zinc Oxide Ointment.....	43
White Petrolatum .....	43
LanaShield Skin Protectant Cream.....	44
D-Cerin Advanced Moisturizing Cream .....	45
Quad Wash Cream with Dimethicone.....	45

### HYGIENE


All Purpose Shampoo and Body Wash.....	49
Moisturizing Shampoo and Body Wash .....	49
Perineum Wash.....	50
Trio Bliss Cleansers.....	51
No Rinse Shampoo Caps.....	51
Hydrea Body Lotion.....	52
DynaSilk with Aloe .....	53
Comfort Bath Oil .....	53
Fresh Mist Odor Eliminator.....	55
Adhesive Tape Remover Pad.....	55

DÜRMA+ 



A close-up photograph showing a medical professional wearing blue nitrile gloves applying a white, porous, advanced wound care dressing to a patient's arm. The dressing is being placed over a wound that has some dried blood. The patient's arm is resting on a light-colored surface.

# ADVANCED WOUND CARE




**Polyurethane Backing:** Air-permeable, moisture-resistant, barrier to bacteria and viruses

**Polyurethane Foam Pad:** Draws in exudate, contours to body, suitable for use with compression

## CURAFOAM™ FOAM DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3011</b>	2" × 2"	10	12	A6209
<b>3012</b>	4" × 4¼"	10	12	A6210
<b>3013</b>	6" × 6"	10	12	A6210

- + Maintains a warm, moist and cushioned healing environment
- + Sterile polyurethane foam helps manage exudate in heavily draining wounds
- + Non-adherent design reduces trauma during dressing changes
- + Helps reduce the risk of contamination during wound care
- + Suitable as a primary or secondary dressing for chronic and acute wounds at stages 2, 3, and 4
- + Appropriate for lower extremity ulcers, donor sites, and partial- to full-thickness wounds
- + Available in multiple sizes to meet various clinical needs




**Backing Layer:** Soft, stretchable, non-woven material allows for excellent airflow

**Wound Pad:** Absorbs exudate and protects the wound

**Adhesive Border:** Gentle on skin for minimal irritation

### DYNAFOAM™ WATERPROOF BORDERED FOAM DRESSING



**Backing Layer:** Soft, stretchable, non-woven material allows for excellent airflow

**Wound Pad:** Absorbs exudate and protects the wound

**Adhesive Border:** Gentle on skin for minimal irritation


### DYNALEVIN™ WATERPROOF ADHESIVE BORDERED FOAM DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3037</b>	4" x 4"	10	12	A6212
<b>3038</b>	6" x 6"	10	12	A6213

- + Maintains a warm, moist healing environment to support healing and help reduce pain
- + Multi-layer design includes absorbent foam and a low-adherent contact layer to protect the wound bed
- + Waterproof outer layer helps shield the wound from external fluids and bacterial contaminants
- + Gentle adhesive border supports secure placement with minimal skin irritation
- + High-absorbency dressing used for heavily exuding wounds
- + Primary or secondary dressing for chronic and acute wounds at stages 2, 3, or 4
- + Suitable for lower extremity ulcers, donor sites, and partial- to full-thickness wounds

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3021</b>	4" x 4"	10	12	A6212
<b>3022</b>	6" x 6"	10	12	A6212

- + Made of absorbent polyurethane foam with waterproof outer layer and adhesive border
- + Helps maintain a warm, moist, and comfortable wound environment to support healing and reduce pain
- + Adhesive backing helps keep contaminants out while securing the dressing to skin
- + Medium-absorbency dressing suitable for exuding wounds
- + Can be used as a primary or secondary dressing for stages 2, 3, or 4
- + Suitable for lower extremity ulcers, donor sites, and full-thickness wounds



**Backing Layer:** Soft, stretchable, non-woven material allows for excellent airflow

**Wound Pad:** Absorbs exudate and protects the wound


**Adhesive Border:** Gentle on skin for minimal irritation

## DYNAGUARD™ WATERPROOF COMPOSITE DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3035</b>	4" × 4"	10	12	A6203
<b>3036</b>	6" × 6"	10	12	A6203

- + Multi-layer, medium-absorbency dressing manages moderately exuding wounds while providing a moist healing environment
- + Waterproof outer layer helps protect against external contaminants
- + Inner foam layer absorbs drainage to support a moist environment
- + Low-adherent contact layer protects the wound bed
- + Adhesive non-woven fabric border conforms gently to skin
- + Primary or secondary dressing to manage chronic and acute wounds at stages 2, 3, and 4






## ISLAND DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3045</b>	6" × 3½"	25	6	–
<b>3074</b>	4" × 13¾"	25	6	–
<b>3492</b>	2" × 3½"	50	6	A6254
<b>3493</b>	4" × 4"	25	6	A6254
<b>3494</b>	6" × 6"	25	6	A6254

- + Soft, non-stick dressing with adhesive border for secure placement
- + Provides excellent absorption and cushioning protection
- + Ideal for use as a cover dressing after packing large, deep wounds
- + Individually packed in peel-open pouches for optimal hygiene and convenience
- + Not made with natural rubber latex



## FOAMFLEX™ NON-ADHESIVE WATERPROOF FOAM DRESSINGS


ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3023</b>	2" × 2"	10	12	A6209
<b>3024</b>	4" × 4"	10	12	A6209
<b>3025</b>	6" × 6"	10	12	A6210

- + Helps maintain a warm, moist wound environment to support healing and reduce pain
- + Semi-permeable film is bonded to the foam, allowing oxygen flow while protecting against contaminants
- + Flexible and easy to conform to hard-to-dress areas
- + High-absorbency design suitable for managing heavily exuding wounds
- + Can be used as primary or secondary dressing for chronic and acute wounds at stages 2, 3, and 4
- + Suitable for lower extremity ulcers, donor sites, and partial- to full-thickness wounds

NON-ADHERENT



SELF-ADHERENT




## DYNASORB™ SUPER ABSORBENT DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>NON-ADHERENT</b>				
<b>3087</b>	4" x 4"	10	12	A6196
<b>3088</b>	6" x 6"	10	12	A6197
<b>SELF-ADHERENT</b>				
<b>3089</b>	4" x 4"	10	12	A6196
<b>3090</b>	6" x 6"	10	12	A6197

- + Multi-layered dressing features a contact layer and waterproof backing to help prevent contamination
- + Super absorbent polymer core rapidly draws in exudate and locks it away to reduce the risk of maceration
- + Soft, conformable construction allows easy application to curved or hard-to-dress areas
- + High absorbency dressing used for heavily exuding wounds
- + Suitable for use as primary or secondary dressing for partial- to full-thickness wounds at stages 2, 3 or 4
- + Appropriate for managing leg and pressure ulcers, surgical wounds, and donor sites
- + Pliable design supports range-of-motion during treatment



\*Non-woven wicking layer and silicone contact layer in self-adherent dressings only.





## SILIGENTLE™ SILICONE FOAM DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCP
<b>NON-BORDERED</b>				
3055	2" × 2"	10	12	A6209
3056	4" × 4"	10	12	A6209
3057	6" × 6"	10	12	A6210
<b>BORDERED</b>				
3058	2" × 2"	10	12	A6212
3059	4" × 4"	10	12	A6212
3060	6" × 6"	10	12	A6212
3061	10" × 9" - Heel	5	12	A6213
3062	7" × 7" - Sacral	5	12	A6212
3063	9" × 9" - Sacral	5	12	A6213
3076	4" × 8"	10	12	-
3091	8" × 6" - Elbow	5	12	-

- + Helps maintain a warm, moist wound environment to support healing
- + Vapor-permeable waterproof film helps protect the wound while allowing moisture balance
- + Soft, comfortable dressing with high absorbency to manage exudate
- + Gentle silicone contact layer adheres securely and moves with the skin, helping to minimize trauma during removal
- + Ideal for high-exuding wounds such as leg ulcers, pressure ulcers, donor sites, and burns
- + Can be used as a primary or secondary dressing for wounds at stages 2, 3, or 4
- + Available in a variety of sizes, with or without a silicone adhesive border






## SILIGENTLE™ AG SILICONE FOAM DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>NON-BORDERED</b>				
<b>3080</b>	2" × 2"	10	5	A6209
<b>3082</b>	4" × 5"	10	5	A6210
<b>3083</b>	6" × 6"	10	5	A6210
<b>BORDERED</b>				
<b>3084</b>	2" × 2"	10	5	A6212
<b>3085</b>	4" × 4"	10	5	A6212
<b>3086</b>	6" × 6"	10	5	A6213

- + Soft, conformable dressing helps maintain a moist wound environment to support healing
- + Silver foam core absorbs exudate and releases silver ions for targeted antimicrobial effect
- + High-absorbency design ideal for managing heavily exuding wounds
- + Suitable for wounds with delayed healing, active infection, or risk of reinfection
- + Vapor-permeable waterproof film allows flexibility and movement without loss of adhesion
- + Gentle silicone contact layer reduces the risk of skin irritation
- + Effective as a primary or secondary dressing for wounds at stages 2, 3, or 4
- + Appropriate for leg ulcers, pressure ulcers, second-degree burns, donor sites, and partial- to full-thickness wounds
- + Available in a variety of sizes, with bordered and non-bordered options



**Backing Layer:** Soft, stretchable, non-woven material allows for excellent airflow

**Absorbent Pad Infused With Silver Ions**


**Adhesive Border:** Gentle on skin for minimal irritation

## DYNAFOAM™ AG BORDERED FOAM DRESSING

ITEM	SIZE	QTY/BX	BX/CS	HCPC
3049	4" x 4"	10	10	A4649

- + Sterile foam dressing infused with silver ions for continuous antimicrobial protection
- + Soft, conformable adhesive design supports secure placement and user comfort
- + Provides high absorbency to manage moderate to heavy exudate
- + Helps maintain a moist wound environment to support healing
- + Waterproof, vapor-permeable border offers protection against contamination
- + Designed for chronic and acute wounds at stages 2, 3, and 4
- + Can be used as a primary or secondary dressing for wounds with delayed healing or at risk of infection
- + Ideal for pressure injuries, surgical wounds, and first- or second-degree burns






## CURAFOAM™ AG FOAM DRESSING

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3047</b>	4" × 4"	10	10	A4649

- + Sterile, non-adherent foam infused with silver ions for antimicrobial protection
- + Helps absorb exudate while maintaining a warm, moist wound environment
- + Provides continuous antimicrobial protection to help reduce colonization on the dressing
- + Highly absorbent foam suitable for chronic and acute wounds at stages 2, 3, and 4
- + Designed for use as a primary dressing on wounds at risk of infection or reinfection
- + Ideal for pressure ulcers, surgical wounds, and first- and second-degree burns



## OIL EMULSION DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3067</b>	3" × 3"	50	4	A6222
<b>3068</b>	3" × 8"	24	6	A6223
<b>3070</b>	5" × 9"	12	6	A6223

- + Knitted gauze dressings impregnated with an emulsion blend of petrolatum and sunflower oil
- + Helps maintain a hydrated wound bed while minimizing patient discomfort
- + Allows exudate to flow into a secondary dressing to help reduce the risk of skin maceration
- + Non-occlusive primary dressing for low-exuding wounds at stages 2, 3, or 4
- + Suitable for donor sites, lacerations, skin grafts, and abrasions



## PETROLEUM GAUZE DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3039</b>	3" x 9"	50	4	A6223
<b>3040</b>	3" x 9"	12	10	A6223
<b>3041</b>	3" x 18"	12	6	A6223

- + Conformable, non-adherent dressing designed to minimize sticking to the wound
- + Infused with petroleum to maintain a moist wound environment and promote autolytic healing
- + Sterile gauze helps to minimize contamination for optimal hygiene
- + Primary dressing for low-draining wounds such as donor sites, minor burns, and non-infected wounds
- + Appropriate for use on full-thickness wounds at stages 2, 3, and 4
- + Not made with natural rubber latex or DEHP




## XEROFORM™ PETROLATUM GAUZE DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3042</b>	3" x 18"	12	6	–
<b>3044</b>	1" x 8"	50	4	–
<b>3053</b>	4" x 4"	25	4	A6222
<b>3054</b>	5" x 9"	50	4	A6223

- + Sterile, non-adherent, knitted gauze dressings designed to maintain a moist wound environment
- + Impregnated with a blend of petrolatum, 3% bismuth, mineral oil, and paraffin oil
- + Non-occlusive design allows the wound to breathe, reducing the risk of skin and tissue maceration
- + Conforms easily to both deep and shallow wounds
- + Primary dressing on low-exuding wounds at stages 2, 3, or 4
- + 3% bismuth helps reduce wound odor






## DYNAGINATE™ CALCIUM ALGINATE DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3010</b>	4" × 4"	10	12	–
<b>3026</b>	2" × 2"	10	12	A6196
<b>3027</b>	4¼" × 4¼"	10	12	A6197
<b>3028</b>	4" × 8"	5	12	A6197
<b>3029</b>	2 g × 30 cm (Rope)	5	12	A6199

- + Calcium alginate dressing transforms into a fibrous gel upon contact with wound exudate
- + Soft gel contours to the wound bed and does not adhere to tissue during dressing changes
- + Highly absorbent; ideal for heavily exuding wounds
- + Maintains a moist wound environment to support the healing process
- + Suitable as a primary dressing for partial- and full-thickness wounds at stages 3 or 4




## DYNAGINATE™ AG CALCIUM ALGINATE DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3030</b>	2" × 2"	10	10	A6196
<b>3031</b>	4" × 5"	10	10	A6197
<b>3032</b>	4" × 8"	5	10	A6197
<b>3033</b>	¾" × 12" Rope	5	10	A6199

- + Combination of calcium alginate and silver ions transforms into a fibrous gel upon contact with wound exudate
- + Creates a moist wound environment and provides an antimicrobial barrier
- + High-absorbency design helps manage heavily exuding wounds
- + Suitable for infected wounds or those at risk of infection
- + Primary dressing ideal for the management of partial- and full-thickness wounds at stages 3 or 4





## DYNADERM™ HYDROCOLLOID DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>GEL</b>				
<b>3014</b>	2" x 2" - Thin	20	12	A6234
<b>3015</b>	4" x 4" - Thin	10	12	A6234
<b>3016</b>	6" x 6" - Thin	5	12	A6238
<b>3017</b>	6" x 7" - Sacral	5	12	A6238
<b>3018</b>	4" x 4" - X-Thin	10	12	A6234
<b>3019</b>	6" x 6" - X-Thin	5	12	A6235
<b>FOAM</b>				
<b>3020</b>	4" x 4" - Foam	5	12	A6209



- + Thermally insulated to reliably protect a variety of non-infected wounds from contamination
- + Semi-occlusive formula helps to maintain moisture and promote autolytic debridement for optimal healing
- + Premium water-resistant materials help prevent bandage peeling for secure protection
- + Flexible, conformable dressings with tapered edges provide a comfortable, anatomical fit
- + Engineered for safe, long-lasting wear, depending on wound condition and exudate levels
- + Caution: Avoid using on diabetic foot wounds without consulting physician for assessment
- + Not made with natural rubber latex



## UNNA BOOT BANDAGES

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>3453</b>	3" × 10 yd.	1	12	A6456
<b>3454</b>	4" × 10 yd.	1	12	A6456
<b>WITH CALAMINE</b>				
<b>3455</b>	3" × 10 yd.	1	12	A6456
<b>3456</b>	4" × 10 yd.	1	12	A6456

- + Designed to help manage venous stasis leg ulcers, thrombophlebitis, lymphatic edema, sprains, strains, and dislocations
- + Designed to maintain a moist healing environment and enhance patient comfort
- + Formulated with soothing zinc oxide paste
- + Unique knit construction allows for more uniform distribution of paste
- + Individually foil-wrapped for optimal hygiene and convenience
- + Finished edges help to minimize risk of fraying into the wound site
- + Available in plain or calamine options
- + Rolled on an inner core for easier, more consistent application



## VIEW GUARD TRANSPARENT DRESSINGS

ITEM	SIZE	QTY/BX	BX/CS	HCPC
<b>3642</b>	2 $\frac{3}{8}$ " × 2 $\frac{3}{4}$ "	100	4	A6257
<b>3643</b>	4" × 4 $\frac{3}{4}$ "	50	4	A6258
<b>3644</b>	6" × 8"	10	8	A6258

- + Designed to help protect the skin, IV sites, and surgical sites
- + Transparent design provides improved visualization of wound site
- + Semi-permeable dressings allow air in while keeping moisture vapor and fluids out
- + Conforms to body contours for a secure, comfortable fit
- + Sterile, not made with natural rubber latex
- + Suggested use: Leave on for up to 3 hours





## PACKING STRIPS - PLAIN & WITH IODOFORM

ITEM	SIZE	QTY/BTL	BTLS/CS	HCPC
<b>PLAIN</b>				
<b>3401</b>	1/4" × 5 yd.	1	12	A6407
<b>3402</b>	1/2" × 5 yd.	1	12	A6407
<b>3403</b>	1" × 5 yd.	1	12	A6407
<b>3404</b>	2" × 5 yd.	1	12	A6407
<b>IODIFORM</b>				
<b>3411</b>	1/4" × 5 yd.	1	12	A6266
<b>3412</b>	1/2" × 5 yd.	1	12	A6266
<b>3413</b>	1" × 5 yd.	1	12	A6266
<b>3414</b>	2" × 5 yd.	1	12	A6266

- + Sterile strips crafted from 100% cotton materials
- + Ideal for use as nasal packing, drainage wicks, or packing wounds
- + Specially woven gauze is designed to minimize the risk of raveling, linting, or fraying
- + Available plain or impregnated with an iodoform solution






## DRESSING RETENTION SHEET

ITEM	SIZE	QTY/BX	BXS/CS	HPC
<b>3542</b>	2" x 11 yd.	1	5	A4452
<b>3544</b>	4" x 11 yd.	1	5	A4452
<b>3546</b>	6" x 11 yd.	1	5	A4452

- + Low-sensitivity adhesive tape designed for secure, comfortable dressing retention
- + Soft, flexible, non-woven fabric is ideal for securing appliances or dressings
- + Porous structure is water-resistant and can cover an entire primary wound dressing conveniently in one application
- + Linerless roll means faster application with no extra step for removing the paper release liner
- + Not made with natural rubber latex



## VIEW GUARD TRANSPARENT DRESSING ROLLS

ITEM	SIZE	QTY/CS	HPC
<b>3637</b>	8" x 11 yd.	12	A4450
<b>3638</b>	10" x 11 yd.	12	A4450
<b>3639</b>	2" x 11 yd.	12	A4450
<b>3640</b>	4" x 11 yd.	12	A4450
<b>3641</b>	6" x 11 yd.	12	A4450

- + Premium transparent dressing designed to help cover and protect skin after needlework
- + Crafted from flexible materials to conform to contours during movement without stressing skin
- + Roll can be easily cut down into desired shape or size for different needs
- + Ideal for use around tattoos, IV sites, surgical sites, incisions, and wounds
- + Suggested use: Leave on for up to 3 hours



## DYNACOL BOVINE COLLAGEN POWDER

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>3094</b>	1 g Packs	10	10	A6010

- + 100% bovine Type 1 collagen powder designed for effective wound management
- + Forms a soft, biodegradable gel upon contact with wound exudate to help promote an optimal healing environment
- + Conforms easily to irregular wound surfaces without the need for cutting or resizing
- + Promotes a moist wound environment for optimal healing conditions
- + Individually packaged in sterile 1-gram pouches for single-use applications
- + Ideal for managing partial and full-thickness wounds, including ulcers, burns, and surgical wounds
- + Designed for use in conjunction with a secondary dressing




## DYNACOL COLLAGEN WOUND DRESSINGS

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>3096</b>	2" x 2"	10	10	A6020
<b>3097</b>	4" x 4"	10	10	A6021
<b>3098</b>	6" x 6"	5	10	A6022
<b>3099</b>	7" x 7"	5	10	A6023

- + 100% bovine Type I collagen wound dressings designed for partial and full-thickness wounds
- + Biodegradable gel forms upon contact with wound exudate to maintain a moist healing environment
- + Pliable, absorbent dressing easily conforms to irregular wound surfaces for comprehensive wound coverage
- + Sterile, single-use dressings, individually packaged in sterile pouches for optimal hygiene and convenience
- + Suitable for use on pressure ulcers, diabetic ulcers, venous ulcers, partial thickness burns, and surgical wounds
- + Available in multiple sizes to accommodate various wound types





## DYNACOL AG COLLAGEN WOUND DRESSING

ITEM	SIZE	QTY/BX	BXS/CS
<b>3095</b>	2" x 2"	10	10
<b>3100</b>	4" x 4"	10	10
<b>3101</b>	6" x 6"	10	10
<b>3105</b>	7" x 7"	5	10

- + Bovine Type I collagen wound dressings combined with silver chloride designed for partial and full-thickness wounds
- + Soft gel forms upon contact with wound exudate to help maintain a moist healing environment
- + Collagen matrix supports new tissue growth and natural healing
- + Silver chloride helps protect the dressing from microbial contamination
- + Pliable, absorbent material conforms to irregular wound surfaces for effective coverage
- + Sterile, single-use dressings individually packaged for hygiene and convenience
- + Suitable for use on ulcers, donor sites, and surgical wounds
- + Available in multiple sizes to accommodate various wound types



## DYNAGEL™ HYDROGEL IMPREGNATED GAUZE DRESSING

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>3064</b>	2" x 2"	15	12	A6231
<b>3065</b>	4" x 4"	15	12	A6231

- + Sterile gauze dressing impregnated with hydrogel for moist wound packing
- + Designed for dry to minimally exudating wounds to support natural healing
- + Hydrogel formulation hydrates dry or necrotic tissue and aids autolytic debridement
- + Requires a secondary dressing to help maintain moisture balance
- + Suitable for management of partial and full-thickness wounds, including dermal and pressure ulcers at stages 2, 3, and 4



## TRANSPARENT FILM WINDOW DRESSING

ITEM	SIZE	QTY/BX	BXS/CS	H CPC
<b>3647</b>	2 $\frac{3}{8}$ " x 2 $\frac{3}{4}$ "	100	4	A6257
<b>3648</b>	4" x 4 $\frac{3}{4}$ "	50	4	A6257


- + Transparent film dressing designed for versatile wound and catheter site care
- + Breathable film supports oxygen exchange to support the healing process
- + Hypoallergenic adhesive helps to reduce irritation on sensitive skin
- + Flexible dressing conforms to body contours and patient movement for improved comfort
- + Features a frame-style design for secure, precise placement
- + Stretch-and-release removal minimizes skin trauma during dressing changes
- + Water-resistant barrier helps to protect against bacteria and fluid contaminants
- + Available in multiple sizes for customized applications

## DYNARULE™ DISPOSABLE PAPER RULER

ITEM	DETAILS	QTY/PAD	PADS/CS
<b>3006</b>	15 cm (150 mm)	25	10

- + Designed to provide direct and reliable wound measurement
- + Disposable, single-use design helps promote optimal hygiene practices
- + Clearly marked in both centimeters and millimeters for versatile use
- + Flexible design conforms to body contours for more accurate assessment
- + Features ample writing space for patient identification and documentation
- + Suitable for initial wound assessment and tracking healing over time
- + Pocket-sized, lightweight rulers are easy to carry and store
- + Not made with natural rubber latex





## DYNARULE™ BULLSEYE MEASURING GUIDE

ITEM	QTY/BX	BXS/CS
<b>3007</b>	250	5

- + Disposable, transparent guide for measuring wound length, width, and diameter
- + Includes both a straight-edge ruler and a bullseye grid
- + Ruler is calibrated in inches and centimeters, with measurements up to 4" (10 cm)
- + Flexible plastic construction conforms to body contours
- + Designed to support consistent wound assessment and documentation



## RETRACTABLE TAPE MEASURE

ITEM	QTY/BX	BXS/CS
<b>3008</b>	6	24


- + Clearly marked in both inches and centimeters for accurate readings
- + Conveniently sized to fit easily in pockets, purses, and kits
- + Measures up to 72 inches or 182 centimeters
- + Flexible design allows for easier measurement around body contours
- + Retracts quickly into durable plastic case for added protection
- + Ideal for wound measurement, weight management, and sports training



## DYNAGEL SILVER WOUND GEL (RX ONLY)

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1279</b>	1.5 oz.	1	24	A6248
<b>1281</b>	3 oz.	1	24	A6248

- + Sterile, prescription-only hydrogel infused with antimicrobial silver to reduce the risk of bacterial growth
- + Absorbs exudate and adheres securely for up to three days
- + Designed to maintain a moist wound environment to support optimal healing
- + Non-staining and non-stinging formula for improved patient comfort
- + Ideal for managing first- and second-degree burns, ulcers, and surgical wounds



## DYNAGEL SKIN PROTECTANT WITH GLYCERIN

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1280</b>	3 oz.		24	A6248

- + Formulated to hydrate and protect wounds and surrounding skin
- + Helps maintain a moist wound environment to support comfort and healing
- + Contains aloe and vitamin E to help soothe irritation and relieve dryness
- + Aids in softening tissue to encourage autolytic debridement
- + Ideal for minor wounds, abrasions, and superficial injuries





## L-MESITRAN® SOFT

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>3078</b>	1.75 oz.	1	24	A6250

- + Formulated with 40% medical-grade honey, hypo-allergenic lanolin, and vitamins C and E to promote optimal healing of chronic wounds
- + Moisturizing formula helps to restore wound pH balance
- + Anti-microbial gel supports cleansing and autolytic debridement
- + Commonly used with secondary dressings to help minimize adherence to the wound bed
- + Ideal for use on partial- and full-thickness wounds, pressure injuries, burns, and trauma wounds at stages 1–4



## DYNAWOUND WOUND CLEANSER SPRAY

ITEM	SIZE	QTY/CS	HCPC
<b>3004</b>	4 fl. oz.	12	A6260
<b>3005</b>	8 fl. oz.	24	A6260

- + Premium cleansing spray designed to help minimize the risk of infection in minor wounds
- + Ideal for use as a convenient supplement to saline irrigation
- + Formulated to help remove exudate and debris to help promote healing
- + Adjustable spray nozzle allows for controlled application
- + Non-toxic and pH-balanced formula for gentle wound cleansing



## NEGATIVE PRESSURE WOUND THERAPY PUMP

ITEM	QTY/CS	HCPC
3150	1	E2402

- + Designed to efficiently administer negative pressure wound therapy (NPWT)
- + Provides gentle powered suction for effective removal of exudate, infectious materials, and tissue debris
- + Offers variable intermittent and continuous therapy modes to accommodate various treatment protocols
- + Effective for managing chronic, acute, traumatic, sub-acute, and dehisced wounds, as well as pressure and diabetic ulcers, burns, flaps, and grafts
- + Easy-to-navigate LCD screen enhances device usability and setting adjustments
- + Compact and lightweight, ideal for transitioning between various care settings
- + Includes an integrated troubleshooting guide for addressing common issues such as canister full, vacuum leak, and blockage alerts
- + Sleep-mode feature extends battery life, facilitating longer mobile use without recharging
- + Battery runtime of up to 72 hours, depending on settings





### CARRY CASE FOR NPWT PUMP



### WOUND VAC DISPOSABLE CANISTERS

ITEM	QTY/CS
------	--------


**3151** 1

- + Durable carry bag designed to provide increased mobility for NPWT patients
- + Ideally sized to fit Dürma+ Negative Pressure Wound Therapy Pump and canister
- + Expandable zipper accommodates larger canisters when needed
- + Clear front panel provides easy visibility of the control panel
- + Transparent back window allows for monitoring of canister fluid levels
- + Adjustable carry strap offers customization for shoulder, waist, or wheelchair use
- + Hook-and-loop closure allows for quick access without removing the pump
- + Convenient charging port supports pump recharging within the case
- + Compact and lightweight design ensures hassle-free portability


ITEM	SIZE	QTY/CS	HCPC
------	------	--------	------

**3121** 140 mL 20 A7000  
**3122** 400 mL 20 A7000

- + Secure, disposable canisters designed for collecting and storing wound exudates
- + Hermetic-seal latch helps to prevent air leakage and sustains pressure under varying conditions
- + User-friendly design with one-handed release enables efficient setup
- + Innovative dam design keeps the membrane dry, reducing false alarms for full canisters
- + Graduated markings aid in accurate fluid tracking and measurement
- + Pre-assembled with drainage tube and male luer for easy connectivity



## SUCTION DRAINAGE TUBES FOR NEGATIVE PRESSURE WOUND THERAPY DEVICE



## FOAM DRESSING

ITEM	DETAILS	SIZE	QTY/CS	HCPC
<b>3147</b>	Y Connector	10 FR	20	A9900
<b>3148</b>	Straight Connector	10 FR	20	A9900

- + Designed to reliably connect wound dressings to NPWT systems
- + Made from high-quality PVC materials for durability and performance
- + 60-inch tubing length enhances patient mobility during treatment
- + Multi-tiered pagoda-style joint provides a secure, customizable fit with dressing kits
- + Anti-Ruhr connector minimizes leakage for uninterrupted therapy
- + Robert splint allows temporary drainage flow stoppage without changing the dressing
- + 3147: Y connector supports simultaneous treatment of multiple wound sites
- + 3148: Straight connector offers effective single-site therapy options
- + Not made with natural rubber latex

ITEM	SIZE	QTY/CS	HCPC
<b>3126</b>	Small, 4" × 2"	10	A6550
<b>3127</b>	Medium, 7" × 4.7"	10	A6550
<b>3128</b>	Large, 10" × 5.9"	10	A6550

- + Premium-quality foam dressings designed for use with NPWT pump systems
- + Made from high-tensile strength foam for optimal durability
- + Flexible foam can be cut to contour a variety of wound shapes, allowing for an individualized fit
- + Available in multiple sizes to accommodate a variety of wound types
- + Single-use and sterile for optimal hygiene
- + Not made with natural rubber latex





## TRANSPARENT WOUND DRAPE

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>3129</b>	12" x 12"	5	20	A6257

- + Designed to provide effective wound protection during negative pressure wound therapy
- + Clear film allows visual monitoring of the wound without removing the entire dressing
- + Helps to ensure proper placement of the foam dressing over the wound
- + Polyurethane film provides breathable coverage to help maintain an optimal healing environment
- + Single-use and sterile for optimal hygiene
- + Not made with natural rubber latex



## FOAM DRESSING KITS

ITEM	SIZE	QTY/CS	HCPC
<b>3123</b>	Small, 4" x 2"	10	A6550
<b>3124</b>	Medium, 7" x 4.7"	10	A6550
<b>3125</b>	Large, 10" x 5.9"	10	A6550

- + Comprehensive dressing kits designed for efficient and comfortable negative pressure wound therapy (NPWT)
- + Engineered for use with an NPWT pump
- + Transparent polyurethane film drapes provide breathable coverage with strong adherence
- + High-tensile strength foam is durable, flexible, and can be cut to contour a variety of wound shapes and sizes, allowing for an individualized fit
- + Dome-shaped silicone port pad offers durability and user comfort
- + Single-use, sterilized components for optimal hygiene
- + Available in multiple sizes to accommodate a variety of wound types

**Care**  **We**  
Care™





# CLINICAL SKIN CARE



## CALMARA - EXTERNAL ANALGESIC

ITEM	SIZE	QTY/CS
<b>1200</b>	7.5 oz.	24

- + Formulated to provide temporary relief of itching from minor skin irritations
- + Helps soothe discomfort from dry skin, insect bites, poison ivy, poison oak, and poison sumac
- + Contains menthol and camphor for a gentle, cooling sensation
- + Safe for use on minor cuts and abrasions
- + Steroid-free formula suitable for sensitive skin
- + Appropriate for healthcare, long-term care, and home environments



## ANTIFUNGAL CREAMS

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>2% MICONAZOLE NITRATE CREAM</b>				
<b>1230-MN</b>	5 g Packets	144	2	A6250
<b>1233-MN</b>	4 oz. Tube	6	4	A6250
<b>1% CLOTRIMAZOLE USP CREAM</b>				
<b>1231-C</b>	1 oz. Tube	1	72	A6250
<b>1233-C</b>	4 oz. Tube	6	4	A6250

- + Soothing topical cream formulated to ease skin irritation
- + Helps relieve itching, scaling, and discomfort from fungal skin infections
- + Effective for athlete's foot, jock itch, and similar conditions
- + Long-lasting, non-greasy formula for comfortable application
- + 1231-c, 1233-c: Formulated with 1% clotrimazole USP
- + 1230-mn, 1233-mn: Formulated with 2% miconazole nitrate
- + Compact size allows for easy packing in first-aid kits and medicine cabinets
- + Ideal for a variety of healthcare settings, including long-term care and home care



## VITAMINS A&D OINTMENTS

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1150</b>	5 g Packet	144	6	A6250
<b>1152</b>	1 oz. Tube	1	72	A6250
<b>1155</b>	4 oz. Tube	6	12	A6250
<b>1157</b>	15 oz. Jar		12	A6250

- + Formulated with vitamins A and D to protect and soothe dry, irritated skin
- + Designed to effectively relieve irritation caused by minor cuts, burns, and diaper rash
- + Infused with lanolin and lavender oil for improved hydration and pleasant scent



## DYNASHIELD SKIN PROTECTANT BARRIER CREAM

## HYDROCORTISONE CREAM USP

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1193</b>	5 g Packet	144	2	A6250
<b>1195</b>	4 oz. Tube		24	A6250
<b>1196</b>	16 oz. Jar		12	A6250
<b>WITH DIMETHICONE</b>				
<b>1198</b>	5 g Packet	144	2	A6250
<b>1199</b>	4 oz. Tube		24	A6250

- + Contains zinc oxide to help prevent and relieve dry skin and diaper rash
- + Enriched with aloe vera and vitamins A, D, and E to nourish skin
- + Available in convenient 5 g packets, 4 oz. tubes and 16 oz. jars
- + 1198-1199: Active ingredient dimethicone helps skin retain moisture and reduce flaking and cracking
- + Ideal for long-term care facilities, first aid kits, home care, and more

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1137</b>	0.9 g Packet	144	12	A6250
<b>1137-25</b>	0.9 g Packet	25	72	A6250
<b>1137UB-6</b>	0.9 g Packet	6	100	A6250
<b>1139</b>	1 oz. Tube	1	72	A6250

- + Formulated with 1% hydrocortisone USP to effectively relieve minor skin irritations
- + Designed to help ease itching and discomfort caused by skin allergies, eczema, and insect bites
- + Topical cream designed for external use only
- + Ideal for use in hospitals, clinics, long-term care facilities, first aid, and home care
- + 1137: Single-use packets allow for easy portability, kit-packing, and medical storage
- + 1139: Easy-to-open tube packed in convenient, retail-friendly packaging





## TRIPLE ANTIBIOTIC OINTMENTS

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1180</b>	0.5 g Packet	144	12	A6250
<b>1180UB-10</b>	0.5 g Packet	10	100	A6250
<b>1182</b>	0.9 g Packet	144	12	A6250
<b>1184</b>	0.5 oz. Tube	1	72	A6250
<b>1185</b>	1 oz. Tube	1	72	A6250

- + Designed to help reduce the risk of infection in minor burns and abrasions
- + Triple-acting formula helps minimize bacterial growth and promote skin healing
- + 1184 & 1185: Easy-open tubes in retail-friendly packaging for quick application and distribution
- + 1180 & 1182: Convenient, single-use packets allow for easy portability and kit packing



## BACITRACIN OINTMENT USP

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1161</b>	0.9 g Packet	144	12	A6250
<b>1162</b>	0.5 oz. Tube	1	72	A6250
<b>1163</b>	1 oz. Tube	1	72	A6250

- + Helps reduce the risk of infection in burns and abrasions
- + Topical antibiotic ointment ideal for first aid, long-term care, and home care
- + 1161: Single-use packets allow for easy portability and kit-packing
- + 1162 & 1163: Easy-open tubes packed in convenient, retail-friendly packaging



## AQUASOOTHE HYDRATING OINTMENT

ITEM	SIZE	QTY/BX	BXS/CS
<b>1463</b>	3.75 oz. Tube		24
<b>1464</b>	14 oz. Jar		12
<b>1465</b>	0.9 g Packet	144	12

- + Multipurpose ointment helps relieve dry, irritated, or chapped skin and lips
- + Formulated with 41% petrolatum to preserve skin's natural moisture
- + Helps protect against the drying effects of wind and cold weather
- + Fragrance-free, non-irritating formula helps to protect minor cuts, scrapes, and burns
- + Ideal for use in hospitals, long-term care facilities, at home, and in the workplace
- + 1465: Convenient packets are ideal for travel, kit packing, and tattoo aftercare



## CALASOOTHE SKIN PROTECTANT MOISTURE BARRIER CREAM

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1271</b>	3.5 g Packet	144	2	A6250
<b>1273</b>	2.5 oz. Jar		36	A6250
<b>1275</b>	4 oz. Tube	1	24	A6250

- + Formulated with menthol, calamine, and lanolin to promote healthy skin and minimize dryness
- + Creates a moisture barrier to reduce irritation due to perspiration, urine, and wound drainage
- + Helps to soothe the discomfort associated with minor cuts, burns, and scrapes
- + Ideal for long-term care facilities, first aid kits, home care, and more





## ZINC OXIDE OINTMENTS

ITEM	SIZE	QTY/CS	HCPC
<b>1190</b>	1 oz. Tube	72	A6250
<b>1191</b>	2 oz. Tube	72	A6250
<b>1192</b>	15 oz. Jar	12	A6250

- + Zinc-enriched formula designed to soothe minor irritations such as diaper rash and chapped skin
- + Helps provide relief from symptoms related to poison ivy, poison oak, and poison sumac
- + Designed to form a protective barrier to seal out moisture and reduce friction and chafing
- + Ideal for hospitals, long-term care facilities, first-aid kits, and home use



## WHITE PETROLATUM

ITEM	SIZE	QTY/BX	BXS/CS	HCPC
<b>1141</b>	1 oz. Tube	1	72	A6250
<b>1145</b>	4 oz. Tube	6	12	A6250
<b>1147</b>	15 oz. Jar		12	A6250

- + Formulated with petrolatum to help protect and soothe skin
- + Forms a moisture barrier to help prevent diaper rash
- + Sensitive formula is ideal for soothing dry, irritated skin



## LANASHIELD SKIN PROTECTANT CREAM

ITEM	SIZE	QTY/BX	BXS/CS	QTY/CS	HCPC
<b>1260</b>	5 g Packet	144	2		A6250
<b>1262</b>	2.5 oz. Jar			36	A6250
<b>1263</b>	4 oz. Tube			24	A6250
<b>1265</b>	14 oz. Jar			12	A6250

- + Formulated with 50% lanolin to help soothe and protect irritated skin
- + Helps relieve discomfort from denuded skin and minor burns
- + Helps relieve chafed, chapped, or cracked skin and lips
- + Aids in preventing thickened skin from becoming hard and brittle
- + Non-greasy formula helps seal out wetness





### D-CERIN ADVANCED MOISTURIZING CREAM

ITEM	SIZE	QTY/BX	BXS/CS
<b>1473</b>	3.75 oz. Tube		24
<b>1474</b>	16 oz. Jar		12

- + Fragrance-free formula helps protect skin from the drying effects of wind and cold
- + Designed to help moisturize and relieve extremely dry, chapped, or cracked skin
- + Active ingredient petrolatum (33%) allows skin to breathe for a healthy glow



### QUAD WASH CREAM WITH DIMETHICONE

ITEM	SIZE	QTY/CS
<b>1388</b>	4 oz. Bottle	24

- + Gentle cleansing cream supports reliable daily hygiene
- + Dimethicone barrier helps protect the skin during routine care
- + Nourishing formula contributes to maintaining overall comfort
- + Moisturizing action helps reduce dryness in delicate areas
- + Ideal for perineal cleansing and incontinence-related hygiene
- + Suitable for acute care, long-term care, rehabilitation, home health, and other assisted-care settings

DYNACARE





**HYGIENE**





## ALL PURPOSE SHAMPOO AND BODY WASH

ITEM	SIZE	QTY/CS
<b>1393</b>	8 fl. oz. Bottle	48
<b>1394</b>	1 gal. Bottle	4

- + All-in-one shampoo and body wash designed to help clean skin and hair
- + Mild formula richly lathers and is gentle on dry skin
- + Can also be used as a hand soap for optimal versatility
- + Lightly scented with pleasant fragrance
- + Ideal for hospitals, long-term care facilities, and home care



## MOISTURIZING SHAMPOO AND BODY WASH

ITEM	SIZE	QTY/CS
<b>1386</b>	8 fl. oz. Bottle	48
<b>1386-800</b>	800mL / Bag	12
<b>1386-1000</b>	1000mL / Bag	10
<b>1387</b>	1 gal. Bottle	4

- + Soap-free botanical gel designed to cleanse both hair and body
- + Moisturizing formula is gentle on dry skin
- + Leaves skin and hair feeling soft and silky
- + Richly lathering gel features a pleasant tropical scent
- + Ideal for routine bathing, shampooing, and hand-washing
- + 1386-800, 1386-1000: Convenient refill bags designed for use with dispensers



## PERINEUM WASH

ITEM	SIZE	QTY/CS
<b>4850</b>	8 fl. oz. Bottle	48
<b>4851</b>	1 gal. Bottle	4

- + No-rinse formula provides convenient cleansing without water
- + Alcohol-free cleanser is formulated with aloe to help soothe delicate skin and reduce the risk of irritation
- + Ideal for incontinence care and personal hygiene
- + Features a gentle, fresh scent



## TRIO BLISS CLEANSERS

ITEM	SIZE	QTY/CS
<b>1382</b>	7.5 fl. oz. Cleansing Foam	12
<b>1383</b>	1 gal. Cleanser	4
<b>1384</b>	7.5 fl. oz. Cleanser	48


- + Provides a convenient, all-in-one option for body wash, shampoo, and perineal care
- + Rinse-free formula supports care routines where access to water is limited
- + Allows for comfortable hygiene care for in-bed and reduced-mobility patients
- + Helps enhance caregiver efficiency with one versatile cleansing product
- + 1382: cleansing foam offers fast, gentle application
- + Ideal for acute care, long-term care, rehabilitation, home health, and other assisted-care settings



## NO-RINSE SHAMPOO CAPS

ITEM	QTY/BG	BGS/CS
<b>4889</b>	6	5

- + Pre-moistened design cleanses and conditions without water
- + Single-step application simplifies hair-care routines
- + Soft, non-woven cap supports comfortable in-bed use
- + One-size, unisex design accommodates most users
- + Single-use format supports hygienic practice
- + Compact design provides convenient travel use
- + Ideal for acute care, long-term care, home health, and other assisted-care settings
- + Not made with natural rubber latex



## HYDRĚVA BODY LOTION

ITEM	SIZE	QTY/CS
<b>1186</b>	8 fl. oz.	12
<b>1187</b>	16 fl. oz.	12
<b>1188</b>	24 fl. oz.	12

- + Formulated to help soften and comfort rough, dry skin
- + Applies smoothly to hands, face, legs, elbows, heels, and other routine-care areas
- + Light consistency supports comfortable daily application
- + Fragrance-free to accommodate sensitivities
- + Supports daily moisturization practices in long-term care, rehabilitation, post-acute, and home-care settings



## DYNASILK WITH ALOE

ITEM	SIZE	QTY/CS
<b>1390</b>	4 fl. oz. Bottle	96
<b>1391</b>	8 fl. oz. Bottle	48
<b>1391-800</b>	800 mL / Bag	12
<b>1391-1000</b>	1000 mL / Bag	10
<b>1392</b>	1 gal. Bottle	4

- + Aloe vera gel formulation helps moisturize and soothe dry skin
- + Non-greasy lotion provides clean-feel application for both hands and body
- + Features a light, fresh scent
- + Ideal for use in hospitals and long-term care environments
- + 1391-800, 1391-1000: Convenient refill bags designed for use with dispensers



## COMFORT BATH OIL

ITEM	SIZE	QTY/CS
<b>1381</b>	7.5 fl. oz. Bottle	48

- + Emollient bath oil formulated for use with water
- + Provides gentle cleansing while helping maintain natural skin moisture
- + Offers a soothing, softening effect that supports skin comfort
- + Designed for use in healthcare, long-term care, and home environments
- + Suitable for tub or bath use only





## FRESH MIST ODOR ELIMINATOR

ITEM	SIZE	QTY/CS
<b>4845</b>	7.5 oz.	48

- + Helps reduce unwanted odors in shared care environments
- + Fine mist releases a fresh scent throughout the surrounding area
- + Helps keep indoor areas fresh between routine cleaning cycles
- + Suitable for patient rooms, bathrooms, hallways, and other common areas
- + Ideal for long-term care, rehabilitation, home care, and other assisted-care settings



## ADHESIVE TAPE REMOVER PAD

ITEM	QTY/BX	BXS/CS	HCPC
<b>1505</b>	100	10	A4456

- + Non-irritating wipes designed to help effectively remove adhesive tapes with minimal discomfort
- + Gentle, effective formula leaves skin free of tape residue
- + Allows for efficient removal without sticking to hair or skin for improved patient comfort
- + Ideal for hospitals, long-term care, and home care settings



155 Chestnut Ridge Rd, Suite 2202 Montvale, NJ 07645  
1-888-dynarex | [dynarex.com](http://dynarex.com)

