

dynarex  
*Resp-C<sub>2</sub>*™



# CONTENTS

## HUMIDIFICATION & FILTRATION

Heat Moisture Exchanger Filters .....	4
Hygroscopic Condenser Humidifiers HMEF .....	6
Heat Moisture Exchangers .....	8
Heat Moisture Exchangers, Tracheostomy.....	10
Filters, Tubes and Connectors .....	12

## VENT CARE SUCTION SOLUTIONS

Closed Suction Catheters Tracheal.....	14
Closed Suction Catheters Endotracheal .....	16

## OXYGEN THERAPY

Oxygen Masks, Medium-Concentration .....	18
Oxygen Masks, High-Concentration.....	20
Oxygen Masks, Non-Rebreather & 3-in-1 .....	22
Venturi Masks .....	24
CO <sub>2</sub> Masks & Accessories .....	26
Nasal Oxygen Cannulas .....	
Standard Tip.....	28
Flared Tip .....	30
Cushion Tip .....	32
Super Soft Touch, Cushion Tip.....	34
Cushion Tip with Ultra Soft Ear Foam .....	36
Cushion Tip with Advantage Ear Foam .....	38
Cushion Tip .....	40
EtCO <sub>2</sub> Nasal Sampling Cannulas.....	42
Oxygen Supply Tubing.....	44
Oxygen Concentrator & Accessories.....	46
Oxygen Cylinders .....	48
Oxygen Regulators & Cylinder Wrench .....	50
Oxygen Tank Carts and Holders .....	52
Oxygen TherapyAccessories .....	54
Connectors and Adaptors.....	56
Pulse Oximeters.....	58

## AEROSOL THERAPY

Sterile Saline and Water .....	60
Aerosol Masks .....	62
Tracheostomy Mask .....	64
Aerosol Drainage Bags .....	66
Disposable Water Taps.....	67
Small VolumeNebulizer 6 cc Cups .....	68
Spring Loaded T-Adaptors.....	70
Pediatric Compressor Nebulizers .....	72
Compressor Nebulizers .....	74
50 PSI Compressor and Accessories.....	76

## SUCTION SOLUTIONS

Single Suction Catheters.....	78
-------------------------------	----

Suction Catheter Kits.....	80
Suction Catheter Mini-Kits .....	82
Sleeved Suction Catheter Kits .....	84
Suction Connecting Tubing .....	86
Yankauer Suction Handles.....	88
Disposable Suction Canisters.....	90
Outer Suction Canisters & Liners.....	91
Homecare Suction Unit .....	92
Homecare Suction Disposable Canister Kit .....	93
Suction Accessories.....	94
Suction Units .....	96

## TRACHEOSTOMY CARE

Tracheostomy Care Kits .....	98
Tracheostomy Care Kit Components.....	100
Tracheostomy Sponges.....	102
Comfort-Foam Trach Tube Holders.....	104
Breathable-Foam Trach Tube Holders.....	106
Soft-Foam Trach TubeHolders .....	107

## MANUAL RESUSCITATION

Manual Pulmonary Resuscitators (MPR).....	108
Air Cushion Masks with Valve .....	110

## BRONCHIAL/ASTHMA MANAGEMENT

Pulmonary Function Measurement.....	112
-------------------------------------	-----

## AIRWAY MANAGEMENT

Endotracheal Tubes, Uncuffed.....	114
Endotracheal Tubes, Cuffed .....	116
Endotracheal Tube Introducer (Bougie) .....	118
Intubating Stylettes .....	119
Nasopharyngeal Airways & Kits .....	120
Laryngeal Mask Airways .....	122
Oral Airways (Berman Type) .....	124
Disposable Guedel Airways.....	126
CPR Products .....	128

# ABOUT

Introducing Dynarex Resp-O<sub>2</sub>™, our innovative line of respiratory therapy products. Resp-O<sub>2</sub>™ features a comprehensive range of respiratory solutions developed by experts who understand the challenges providers face while administering proper care. Since our founding in 1967, Dynarex has been one of the industry's leaders in disposable and durable medical equipment. As your partner in the marketplace, we continue to diversify our portfolio across the entire continuum of care because the success of your business matters to us. When it comes to doing business with Dynarex, we choose to work like we're family, because We Care Like Family™.



**dynarex**  
*Resp-C<sub>2</sub>*™





# HEAT MOISTURE EXCHANGER FILTERS



## Heat Moisture Exchanger + B/V Filter (HMEF) 1000

CO<sub>2</sub> PORT

**REORDER NO. 30100**

### ADULT

#### TIDAL VOLUME

150-1000 mL

#### RESISTANCE

at flow 30 L/min 0.8 cm H<sub>2</sub>O

at flow 60 L/min 1.9 cm H<sub>2</sub>O

#### DEAD SPACE

70 mL

#### MOISTURE OUTPUT

250 mL 36.0 mg H<sub>2</sub>O/L

500 mL 34.0 mg H<sub>2</sub>O/L

1000 mL 33.0 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.999%

Viral: >99.99%

#### WEIGHT

24.5 g

#### CONNECTIONS

22M/15F-22F/15M

#### PACKAGE QUANTITY

50/Case

Dynarex Resp-O<sub>2</sub>™ HMEFs are single-use respiratory devices designed to provide humidification and help prevent mucosal dryness during mechanical patient ventilation and other respiratory treatments. The integrated bacterial/viral filter reduces the risk of contamination and provides protection from airborne particles and microorganisms.

- Features a low-resistance design to best support patient's breathing
- Lightweight design for maximum patient comfort
- Transparent housing enables effortless placement
- Sampling port proximal to patient adds to overall convenience
- Not made with natural rubber latex or DEHP



### Heat Moisture Exchanger + B/V Filter (HMEF) 1000

CO<sub>2</sub> PORT • FLEX TUBE

**REORDER NO. 30101**

#### ADULT

#### TIDAL VOLUME

150-1000 mL

#### RESISTANCE

at flow 30 L/min 0.8 cm H<sub>2</sub>O  
at flow 60 L/min 1.9 cm H<sub>2</sub>O

#### DEAD SPACE

105 (70) mL

#### MOISTURE OUTPUT

250 mL 36.0 mg H<sub>2</sub>O/L  
500 mL 34.0 mg H<sub>2</sub>O/L  
1000 mL 33.0 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.999%  
Viral: >99.99%

#### WEIGHT

33 g

#### CONNECTIONS

15F-22F/15M

#### PACKAGE QUANTITY

30/Case



### Heat Moisture Exchanger + B/V Filter (HMEF) 1500

CO<sub>2</sub> PORT

**REORDER NO. 30102**

#### ADULT

#### TIDAL VOLUME

150-1500 mL

#### RESISTANCE

at flow 30 L/min 0.95 cm H<sub>2</sub>O  
at flow 60 L/min 1.7 cm H<sub>2</sub>O

#### DEAD SPACE

55 mL

#### MOISTURE OUTPUT

250 mL 39.8 mg H<sub>2</sub>O/L  
500 mL 38.7 mg H<sub>2</sub>O/L  
1000 mL 39.0 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.998%  
Viral: >99.96%

#### WEIGHT

28.5 g

#### CONNECTIONS

22M-15F/22F-15M

#### PACKAGE QUANTITY

50/Case



### Heat Moisture Exchanger + B/V Filter (HMEF) 1000

CO<sub>2</sub> PORT • FLEX TUBE • LOW PROFILE

**REORDER NO. 30103**

#### ADULT

#### TIDAL VOLUME

150-1000 mL

#### RESISTANCE

at flow 30 L/min 0.8 cm H<sub>2</sub>O  
at flow 60 L/min 1.9 cm H<sub>2</sub>O

#### DEAD SPACE

83 (43) mL

#### MOISTURE OUTPUT

250 mL 35.0 mg H<sub>2</sub>O/L  
500 mL 33.0 mg H<sub>2</sub>O/L  
1000 mL 31.0 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.999%  
Viral: >99.99%

#### WEIGHT

30.5 g

#### CONNECTIONS

15F-22F/15M

#### PACKAGE QUANTITY

30/Case



### Heat Moisture Exchanger + B/V Filter (HMEF) 1000

CO<sub>2</sub> PORT • LOW PROFILE

**REORDER NO. 30104**

#### ADULT

#### TIDAL VOLUME

150-1000 mL

#### RESISTANCE

at flow 30 L/min 0.8 cm H<sub>2</sub>O  
at flow 60 L/min 1.9 cm H<sub>2</sub>O

#### DEAD SPACE

43 mL

#### MOISTURE OUTPUT

250 mL 35.0 mg H<sub>2</sub>O/L  
500 mL 33.0 mg H<sub>2</sub>O/L  
1000 mL 31.0 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.999%  
Viral: >99.99%

#### WEIGHT

21 g

#### CONNECTIONS

22M/15F-22F/15M

#### PACKAGE QUANTITY

50/Case



### Heat Moisture Exchanger (HMEF) 1500 + (HEPA) Filter

CO<sub>2</sub> PORT

**REORDER NO. 30112**

#### ADULT

#### TIDAL VOLUME

250-1500 mL

#### RESISTANCE

at flow 30 L/min 1.5 cm H<sub>2</sub>O  
at flow 60 L/min 2.8 cm H<sub>2</sub>O

#### DEAD SPACE

85 mL

#### MOISTURE OUTPUT

250 mL 37.3 mg H<sub>2</sub>O/L  
500 mL 36.0 mg H<sub>2</sub>O/L  
1000 mL 32.3 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.9999%  
Viral: >99.999%

#### WEIGHT

41 g

#### CONNECTIONS

22M/15F-15M

#### PACKAGE QUANTITY

25/Case



# HYGROSCOPIC CONDENSER HUMIDIFIERS HMEF



## Hygrosopic Condenser Humidifier HMEF 1500

B/V FILTER • CO<sub>2</sub> PORT

**REORDER NO. 30105**

### ADULT

#### TIDAL VOLUME

250-1500 mL

#### RESISTANCE

at flow 30 L/min 0.8 cm H<sub>2</sub>O

at flow 60 L/min 1.8 cm H<sub>2</sub>O

#### DEAD SPACE

60 mL

#### MOISTURE OUTPUT

250 mL 37.3 mg H<sub>2</sub>O/L

500 mL 33.3 mg H<sub>2</sub>O/L

1000 mL 31.3 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.999%

Viral: >99.99%

#### WEIGHT

36 g

#### CONNECTIONS

22M/15F-22F/15M

#### PACKAGE QUANTITY

25/Case



## Hygrosopic Condenser Humidifier HMEF 1500

B/V FILTER • CO<sub>2</sub> PORT • FLEX TUBE

**REORDER NO. 30106**

### ADULT

#### TIDAL VOLUME

250-1500 mL

#### RESISTANCE

at flow 30 L/min 0.9 cm H<sub>2</sub>O

at flow 60 L/min 2.0 cm H<sub>2</sub>O

#### DEAD SPACE

89 (60) mL

#### MOISTURE OUTPUT

250 mL 37.3 mg H<sub>2</sub>O/L

500 mL 33.3 mg H<sub>2</sub>O/L

1000 mL 31.3 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.999%

Viral: >99.99%

#### WEIGHT

45 g

#### CONNECTIONS

15F-22F/15M

#### PACKAGE QUANTITY

25/Case



## Hygrosopic Condenser Humidifier HMEF 600

B/V FILTER • CO<sub>2</sub> PORT

**REORDER NO. 30107**

### PEDIATRIC

#### TIDAL VOLUME

70-600 mL

#### RESISTANCE

at flow 15 L/min 0.8 cm H<sub>2</sub>O

at flow 30 L/min 1.6 cm H<sub>2</sub>O

#### DEAD SPACE

19 mL

#### MOISTURE OUTPUT

250 mL 33.3 mg H<sub>2</sub>O/L

500 mL 27.8 mg H<sub>2</sub>O/L

#### FILTRATION EFFICIENCY

Bacterial: >99.999%

Viral: >99.99%

#### WEIGHT

19 g

#### CONNECTIONS

22M/15F-15M

#### PACKAGE QUANTITY

25/Case

Dynarex Resp-O<sub>2</sub>™ Hygrosopic HMEFs are single-use respiratory devices designed to provide humidification and help prevent mucosal dryness during mechanical patient ventilation and other respiratory treatments. The integrated bacterial/viral filter reduces the risk of contamination and provides protection from airborne particles and microorganisms.

- Features a low-resistance design to best support patient's breathing
- Lightweight design for maximum patient comfort
- Transparent housing enables effortless placement
- Sampling port proximal to patient adds to overall convenience
- Not made with natural rubber latex or DEHP





# HEAT MOISTURE EXCHANGERS



Dynarex Resp-O<sub>2</sub>™ HMEs are single-use respiratory devices designed to provide humidification and help prevent mucosal dryness during mechanical patient ventilation and other respiratory treatments.

- Features a low-resistance design to best support patient's breathing
- Lightweight design for maximum patient comfort
- Transparent housing enables effortless placement
- 30203-30205: Sampling port proximal to patient adds to overall convenience
- Not made with natural rubber latex or DEHP



### Heat Moisture Exchanger (HME) 1500

**REORDER NO. 30201**

#### ADULT

#### TIDAL VOLUME

250-1500 mL

#### RESISTANCE

at flow 30 L/min 1.0 cm H<sub>2</sub>O  
at flow 60 L/min 2.5 cm H<sub>2</sub>O

#### DEAD SPACE

38 mL

#### MOISTURE OUTPUT

1000 mL 30.0 mg H<sub>2</sub>O/L

#### WEIGHT

9 g

#### CONNECTIONS

Patient Side: 22M/15F  
Machine Side: 15M

#### PACKAGE QUANTITY

50/Case



### Heat Moisture Exchanger (HME) 1500

INTEGRATED FLEX TUBE

**REORDER NO. 30202**

#### ADULT/PEDIATRIC

#### TIDAL VOLUME

250-1500 mL

#### RESISTANCE

at flow 30 L/min 0.5 cm H<sub>2</sub>O  
at flow 60 L/min 1.4 cm H<sub>2</sub>O

#### DEAD SPACE

90 mL

#### MOISTURE OUTPUT

1000 mL 30.0 mg H<sub>2</sub>O/L

#### WEIGHT

15.5 g

#### CONNECTIONS

15F-22F

#### PACKAGE QUANTITY

50/Case



### Heat Moisture Exchanger (HME) 1500

HYGROSCOPIC CONDENSER  
HUMIDIFIER • CO<sub>2</sub> PORT

**REORDER NO. 30203**

#### ADULT

#### TIDAL VOLUME

150-1500 mL

#### RESISTANCE

at flow 30 L/min 1.0 cm H<sub>2</sub>O  
at flow 60 L/min 2.8 cm H<sub>2</sub>O

#### DEAD SPACE

34 mL

#### MOISTURE OUTPUT

250 mL 37.7 mg H<sub>2</sub>O/L  
500 mL 35.5 mg H<sub>2</sub>O/L  
1000 mL 33.5 mg H<sub>2</sub>O/L

#### WEIGHT

15 g

#### CONNECTIONS

22M/15F-15M

#### PACKAGE QUANTITY

50/Case



### Heat Moisture Exchanger (HME) 1500

HYGROSCOPIC CONDENSER  
HUMIDIFIER • CO<sub>2</sub> PORT • FLEX TUBE

**REORDER NO. 30204**

#### ADULT

#### TIDAL VOLUME

150-1500 mL

#### RESISTANCE

at flow 30 L/min 1.1 cm H<sub>2</sub>O  
at flow 60 L/min 3.0 cm H<sub>2</sub>O

#### DEAD SPACE

62 (34) mL

#### MOISTURE OUTPUT

250 mL 37.7 mg H<sub>2</sub>O/L  
500 mL 35.7 mg H<sub>2</sub>O/L  
1000 mL 33.5 mg H<sub>2</sub>O/L

#### WEIGHT

24 g

#### CONNECTIONS

15F-15M

#### PACKAGE QUANTITY

30/Case



### Heat Moisture Exchanger (HME) 600

HYGROSCOPIC CONDENSER  
HUMIDIFIER • CO<sub>2</sub> PORT

**REORDER NO. 30205**

#### PEDIATRIC

#### TIDAL VOLUME

70-600 mL

#### RESISTANCE

at flow 30 L/min 0.8 cm H<sub>2</sub>O  
at flow 60 L/min 2.2 cm H<sub>2</sub>O

#### DEAD SPACE

21 mL

#### MOISTURE OUTPUT

250 mL 34.5 mg H<sub>2</sub>O/L  
500 mL 32.5 mg H<sub>2</sub>O/L

#### WEIGHT

9 g

#### CONNECTIONS

22M/15F-15M

#### PACKAGE QUANTITY

50/Case

# HEAT MOISTURE EXCHANGERS

## TRACHEOSTOMY



Heat Moisture Exchanger  
(HME) 1000 Tracheostomy

INTEGRAL OXYGEN PORT

REORDER NO. 30300

ADULT

TIDAL VOLUME

150-1250 mL

RESISTANCE

at flow 60 L/min 0.2 cm H<sub>2</sub>O

DEAD SPACE

16 mL

MOISTURE OUTPUT

500 mL 30.0 mg H<sub>2</sub>O/L

WEIGHT

5.6 g

CONNECTIONS

15F

PACKAGE QUANTITY

50/Case

10



Heat Moisture Exchanger  
(HME) 1250 Tracheostomy

T-SHAPE • SUCTION PORT

REORDER NO. 30301

ADULT

TIDAL VOLUME

150-1250 mL

RESISTANCE

at flow 30 L/min 0.2 cm H<sub>2</sub>O

DEAD SPACE

15 mL

MOISTURE OUTPUT

250 mL 36.3 mg H<sub>2</sub>O/L

500 mL 30.0 mg H<sub>2</sub>O/L

WEIGHT

5.5 g

CONNECTIONS

15F

PACKAGE QUANTITY

50/Case



Heat Moisture Exchanger  
(HME) 1250 Tracheostomy

T-SHAPE • SUCTION PORT  
FRONT OXYGEN PORT

REORDER NO. 30302

ADULT

TIDAL VOLUME

150-1250 mL

RESISTANCE

at flow 30 L/min 0.2 cm H<sub>2</sub>O

DEAD SPACE

17 mL

MOISTURE OUTPUT

250 mL 36.6 mg H<sub>2</sub>O/L

500 mL 33.4 mg H<sub>2</sub>O/L

WEIGHT

4.5 g

CONNECTIONS

15F

PACKAGE QUANTITY

50/Case



Heat Moisture Exchanger  
(HME) 1500 Tracheostomy

T-SHAPE • SWIVEL OXYGEN PORT

REORDER NO. 30303

ADULT

TIDAL VOLUME

75-1500 mL

RESISTANCE

at flow 30 L/min 0.15 cm H<sub>2</sub>O

DEAD SPACE

15 mL

MOISTURE OUTPUT

250 mL 36.8 mg H<sub>2</sub>O/L

500 mL 36.2 mg H<sub>2</sub>O/L

WEIGHT

7.8 g

CONNECTIONS

15F

PACKAGE QUANTITY

50/Case



Heat Moisture Exchanger  
(HME) 1500 Tracheostomy

BOW TIE

REORDER NO. 30305

ADULT

TIDAL VOLUME

75-1500 mL

RESISTANCE

at flow 30 L/min 0.2 cm H<sub>2</sub>O

DEAD SPACE

12 mL

MOISTURE OUTPUT

250 mL 34.9 mg H<sub>2</sub>O/L

500 mL 32.9 mg H<sub>2</sub>O/L

WEIGHT

4.8 g

CONNECTIONS

15F

PACKAGE QUANTITY

50/Case



Dynarex Resp-O<sub>2</sub>™ Trach HMEs are single-use, disposable respiratory devices designed to provide humidification for tracheostomy patients.

- Features a low-resistance design to best support patient's breathing
- Lightweight design for maximum patient comfort
- Transparent housing enables effortless placement
- Not made with natural rubber latex or DEHP



### Heat Moisture Exchanger (HME) Tracheostomy

SIDE OXYGEN PORT

**REORDER NO. 30309**

#### ADULT/PEDIATRIC

##### TIDAL VOLUME

300-1500 mL

##### RESISTANCE

at flow 30 L/min 0.8 cm H<sub>2</sub>O  
at flow 60 L/min 1.8 cm H<sub>2</sub>O

##### DEAD SPACE

15 mL

##### MOISTURE OUTPUT

250 mL 32.7 mg H<sub>2</sub>O/L  
500 mL 28.6 mg H<sub>2</sub>O/L  
1000 mL 25.0 mg H<sub>2</sub>O/L

##### WEIGHT

8.4 g

##### CONNECTIONS

15F

##### PACKAGE QUANTITY

50/Case



### Pediatric/Neonatal Mini Tracheostomy (HME) 50

**REORDER NO. 30210**

#### PEDIATRIC/NEONATAL

##### TIDAL VOLUME

5-50 mL

##### RESISTANCE

at flow 15 L/min 0.4 cm H<sub>2</sub>O  
at flow 30 L/min 2.0 cm H<sub>2</sub>O

##### DEAD SPACE

2.2 mL

##### MOISTURE OUTPUT

50 mL 14.4 mg H<sub>2</sub>O/L

##### WEIGHT

1.6 g

##### CONNECTIONS

15F-15M

##### PACKAGE QUANTITY

50/Case



**dynarex**  
**Resp-O<sub>2</sub>**



# FILTERS, TUBES AND CONNECTORS



## Advantage Filter

UNIVERSAL RESPIRATORY  
B/V FILTER

**REORDER NO. 30400**

### TIDAL VOLUME

150-1500 mL

### RESISTANCE

at flow 30 L/min 0.68 cm H<sub>2</sub>O  
at flow 60 L/min 1.4 cm H<sub>2</sub>O

### DEAD SPACE

45 mL

### FILTRATION EFFICIENCY

Bacterial: >99.998%  
Viral: >99.96%

### WEIGHT

21 g

### CONNECTIONS

Patient Side: 22M/15F  
Machine Side: 22F/15M

### PACKAGE QUANTITY

50/Case



## (HEPA) Filter

CO<sub>2</sub> PORT

**REORDER NO. 30401**

### TIDAL VOLUME

250-1500 mL

### RESISTANCE

at flow 30 L/min 1.2 cm H<sub>2</sub>O  
at flow 60 L/min 2.6 cm H<sub>2</sub>O

### DEAD SPACE

68 mL

### FILTRATION EFFICIENCY

Bacterial: >99.9999%  
Viral: >99.999%

### WEIGHT

39 g

### CONNECTIONS

22M/15F-22F/15M

### PACKAGE QUANTITY

25/Case

## Filters

Dynarex Resp-O<sub>2</sub>™ Bacterial/Viral Filters are designed to improve the safety of your facility and patients on mechanical ventilation by reducing the risk of inhaling harmful pathogens.

- Efficiently protects patient and equipment from cross-contamination
- Not made with natural rubber latex or DEHP



## Universal 6" Flex Tube

**REORDER NO. 30600**

### ADULT/PEDIATRIC

### DEAD SPACE

40 mL

### WEIGHT

9.5 g

### CONNECTIONS

15F/22F

### PACKAGE QUANTITY

50/Case



## Tracheostomy Adult Flex Connector

15 MM OD × 22 MM OD/15 MM ID

**REORDER NO. 30601**

### ADULT

### DEAD SPACE

8 mL

### WEIGHT

6.5 g

### CONNECTIONS

22M/15F-15M

### PACKAGE QUANTITY

50/Case



## Tracheostomy Pediatric Flex Connector

15 MM OD × 15 MM ID

**REORDER NO. 30602**

### PEDIATRIC

### DEAD SPACE

8 mL

### WEIGHT

5.5 g

### CONNECTIONS

15F-15M

### PACKAGE QUANTITY

50/Case

## Tubes/Connectors

Dynarex Resp-O<sub>2</sub>™ Flex Tube/Flex Connectors add length to vertical circuits and HMEs to help improve system flexibility and range of motion.

- Provides added flexibility to ventilator circuits
- Lightweight, pliable design
- Not made with natural rubber latex or DEHP








# CLOSED SUCTION CATHETERS

## TRACHEAL

Length: 30.5 cm/12 in



T-PIECE

Transparent T-Piece facilitates visualization of the pathway

DOUBLE SWIVEL ELBOW

promotes greater flexibility and movement

MDI Port

One-way Irrigation Port helps maintain ventilator circuit pressure and accepts saline injections for cleaning

10 Fr

T-Piece

REORDER NO. 31000

#### PEDIATRIC

#### ACCESSORIES

MDI Adaptor  
15 mm Swivel Adaptor  
22 mm Flex Tube  
15 mm × 22 mm Adaptor  
Change of Day

PACKAGE QUANTITY  
20/Case

12 Fr

T-Piece

REORDER NO. 31001

#### ADULT/PEDIATRIC

#### ACCESSORIES

MDI Adaptor  
15 mm Swivel Adaptor  
22 mm Flex Tube  
15 mm × 22 mm Adaptor  
Change of Day

PACKAGE QUANTITY  
20/Case

14 Fr

T-Piece

REORDER NO. 31002

#### ADULT

#### ACCESSORIES

MDI Adaptor  
15 mm Swivel Adaptor  
22 mm Flex Tube  
15 mm × 22 mm Adaptor  
Change of Day

PACKAGE QUANTITY  
20/Case

16 Fr

T-Piece

REORDER NO. 31003

#### ADULT

#### ACCESSORIES

MDI Adaptor  
15 mm Swivel Adaptor  
22 mm Flex Tube  
15 mm × 22 mm Adaptor  
Change of Day

PACKAGE QUANTITY  
20/Case

**dynarex**  
**Resp-Q<sub>2</sub>**





Dynarex Resp-O<sub>2</sub>™ Closed Suction Trach Catheters are disposable respiratory devices designed to safely remove secretions from the respiratory tract of ventilator-dependent patients. The device allows ventilation to continue during suctioning, helps prevent cross-contamination, and protects caregivers from airborne pathogens.

- Protects patient from contact with outside pathogens and prevents contamination of the surrounding environment
- Smooth, flexible tubing minimizes patient discomfort
- Pressure wiper seal helps maintain PEEP
- Flex Tube provides adjustment and flexibility to the ventilator circuit
- Designed for 24-hr, single-patient use
- Not made with DEHP

**10 Fr**

Double Swivel Elbow

**REORDER NO. 31010**

**PEDIATRIC**

**ACCESSORIES**

MDI Port  
15 mm × 22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**12 Fr**

Double Swivel Elbow

**REORDER NO. 31011**

**ADULT/PEDIATRIC**

**ACCESSORIES**

MDI Port  
15 mm × 22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**14 Fr**

Double Swivel Elbow

**REORDER NO. 31012**

**ADULT**

**ACCESSORIES**

MDI Port  
15 mm × 22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**16 Fr**

Double Swivel Elbow

**REORDER NO. 31013**

**ADULT**

**ACCESSORIES**

MDI Port  
15 mm × 22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case



**DOUBLE SWIVEL ELBOW ACCESSORIES**




**T-PIECE ACCESSORIES**



# CLOSED SUCTION CATHETERS ENDOTRACHEAL

Length: 54.0 cm/21.3 in



**10 Fr**

T-Piece

**REORDER NO. 31020**

**PEDIATRIC**

**ACCESSORIES**

22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**12 Fr**

T-Piece

**REORDER NO. 31021**

**ADULT/PEDIATRIC**

**ACCESSORIES**

22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**14 Fr**

T-Piece

**REORDER NO. 31022**

**ADULT**

**ACCESSORIES**

22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**16 Fr**

T-Piece

**REORDER NO. 31023**


**ADULT**

**ACCESSORIES**

22 mm Flex Tube  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**dynarex**  
**Resp-C<sub>2</sub>**



Dynarex Resp-O<sub>2</sub>™ Closed Suction Endotracheal Catheters are disposable respiratory devices designed to safely remove secretions from the respiratory tract of ventilator-dependent patients. The device allows ventilation to continue during suctioning, helps prevent cross-contamination, and protects caregivers from airborne pathogens.

- Protects patient from contact with outside pathogens and prevents contamination of the surrounding environment
- Smooth, flexible tubing minimizes patient discomfort
- Pressure wiper seal helps maintain PEEP
- Flex Tube provides adjustment and flexibility to the ventilator circuit
- Designed for 24-hr, single-patient use
- Not made with DEHP



**10 Fr**

Double Swivel Elbow

**REORDER NO. 31030**

**PEDIATRIC**

**ACCESSORIES**

15 mm x 22 mm Flex Tube  
MDI Port  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**12 Fr**

Double Swivel Elbow

**REORDER NO. 31031**

**ADULT/PEDIATRIC**

**ACCESSORIES**

15 mm x 22 mm Flex Tube  
MDI Port  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**14 Fr**

Double Swivel Elbow

**REORDER NO. 31032**

**ADULT**

**ACCESSORIES**

15 mm x 22 mm Flex Tube  
MDI Port  
Change of Day

**PACKAGE QUANTITY**  
20/Case

**16 Fr**

Double Swivel Elbow

**REORDER NO. 31033**

**ADULT**

**ACCESSORIES**

15 mm x 22 mm Flex Tube  
MDI Port  
Change of Day

**PACKAGE QUANTITY**  
20/Case



# OXYGEN MASKS

Medium-Concentration



STANDARD



INFANT MASK



ELONGATED

Dynarex Resp-O<sub>2</sub>™ Medium-Concentration Masks are designed to effectively deliver oxygen at a concentration rate of 40% to 50%. These masks come with light green oxygen supply tubing, designed to help improve tube visibility.

- Transparent vinyl mask features smooth edges
- Elastic strap allows for easy adjustment
- Adult and pediatric masks feature an adjustable nose clip
- Includes 7' of crush-resistant, star lumen oxygen supply tubing
- Not made with natural rubber latex or DEHP



Includes  
7' (2.1 m)  
Crush  
Resistant  
Tubing

### Elongated Mask

STANDARD CONNECTOR

**REORDER NO. 33005**

ADULT

PACKAGE QUANTITY  
50/Case

### Elongated Mask

UNIVERSAL CONNECTOR

**REORDER NO. 33006**

ADULT

PACKAGE QUANTITY  
50/Case

### Elongated Mask

STANDARD CONNECTOR

**REORDER NO. 33007**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### Elongated Mask

UNIVERSAL CONNECTOR

**REORDER NO. 33008**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### Standard Mask

STANDARD CONNECTOR

**REORDER NO. 33000**

ADULT

PACKAGE QUANTITY  
50/Case

### Standard Mask

UNIVERSAL CONNECTOR

**REORDER NO. 33001**

ADULT

PACKAGE QUANTITY  
50/Case

### Standard Mask

STANDARD CONNECTOR

**REORDER NO. 33002**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### Standard Mask

UNIVERSAL CONNECTOR

**REORDER NO. 33003**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### Standard Mask

STANDARD CONNECTOR

**REORDER NO. 33004**

INFANT

PACKAGE QUANTITY  
50/Case



# OXYGEN MASKS

## High-Concentration



### Standard Mask

1000 ML RESERVOIR BAG  
STANDARD CONNECTOR

REORDER NO. 33009

ADULT

PACKAGE QUANTITY  
50/Case

### Standard Mask

600 ML RESERVOIR BAG  
STANDARD CONNECTOR

REORDER NO. 33010

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### Standard Mask

600 ML RESERVOIR BAG  
STANDARD CONNECTOR

REORDER NO. 33011

INFANT

PACKAGE QUANTITY  
50/Case

### Elongated Mask

1000 ML RESERVOIR BAG  
STANDARD CONNECTOR

REORDER NO. 33012

ADULT

PACKAGE QUANTITY  
50/Case

### Elongated Mask

600 ML RESERVOIR BAG  
STANDARD CONNECTOR

REORDER NO. 33013

PEDIATRIC

PACKAGE QUANTITY  
50/Case

Star Lumen

Includes  
7' (2.1 m)  
Crush  
Resistant  
Tubing

Dynarex Resp-O<sub>2</sub>™ High-Concentration Masks are designed to effectively deliver oxygen at concentrations of 80%-100%. The masks feature a reservoir bag to ensure sufficient oxygen supply and light green oxygen supply tubing, designed to help improve tube visibility.

- Transparent vinyl mask features smooth edges
- Mask features an adjustable nose clip and elastic strap
- Includes a vinyl reservoir bag and 7' of crush-resistant, star lumen oxygen supply tubing
- Not made with natural rubber latex or DEHP



**dynarex**  
**Resp-O<sub>2</sub>**



# OXYGEN MASKS

## Non-Rebreather & 3-in-1

Dynarex Resp-O<sub>2</sub>™ Non-Rebreather Masks are designed to effectively deliver oxygen at high concentration rates in emergency situations. The masks feature a safety valve and a reservoir bag to ensure sufficient oxygen supply. Dynarex also offers a 3-in-1 mask that can be used as a medium-concentration, high-concentration, or non-rebreathing mask.

Includes  
7' (2.1 m)  
Crush  
Resistant  
Tubing

- Transparent vinyl mask features smooth edges
- Mask features an adjustable nose clip and elastic strap
- Includes a vinyl reservoir bag and 7' of crush-resistant, star lumen oxygen supply tubing
- Not made with natural rubber latex or DEHP

**dynarex**  
**Resp-O<sub>2</sub>**


**33055: 3-IN-1**

### 3-in-1 Elongated Mask

1000 ML RESERVOIR BAG  
STANDARD CONNECTOR

**REORDER NO. 33055**

ADULT

PACKAGE QUANTITY  
50/Case


**33054: INFANT MASK**


Star Lumen

### Non-Rebreather Elongated Mask with Safety Vent

1000 ML RESERVOIR BAG  
STANDARD CONNECTOR

**REORDER NO. 33050**

ADULT

PACKAGE QUANTITY  
50/Case

### Non-Rebreather Elongated Mask with Safety Vent

600 ML RESERVOIR BAG  
STANDARD CONNECTOR

**REORDER NO. 33051**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### Non-Rebreather Elongated Mask

1000 ML RESERVOIR BAG  
STANDARD CONNECTOR

**REORDER NO. 33052**

ADULT

PACKAGE QUANTITY  
50/Case

### Non-Rebreather Elongated Mask

600 ML RESERVOIR BAG  
STANDARD CONNECTOR

**REORDER NO. 33053**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### Non-Rebreather Standard Mask

600 ML RESERVOIR BAG  
STANDARD CONNECTOR

**REORDER NO. 33054**

INFANT

PACKAGE QUANTITY  
50/Case




# VENTURI MASKS

Dynarex Resp-O<sub>2</sub>™ Venturi Masks are designed to deliver oxygen at fixed concentrations ranging from 24% to 50%. These masks come with light green oxygen supply tubing, designed to help improve tube visibility.

- Transparent, molded vinyl mask features smooth edges
- Mask features an adjustable nose clip and elastic strap
- Includes 7' of crush-resistant, star lumen oxygen supply tubing
- Not made with natural rubber latex or DEHP





33070 & 33071 DILUTORS



33072 & 33073 ADJUSTABLE DIAL



Includes  
7' (2.1 m)  
Crush  
Resistant  
Tubing

Select A Vent  
Elongated Mask  
with 6 Dilutors

STANDARD CONNECTOR

REORDER NO. 33070

ADULT  
.....  
PACKAGE QUANTITY  
50/Case

Select A Vent  
Elongated Mask  
with 6 Dilutors

STANDARD CONNECTOR

REORDER NO. 33071

PEDIATRIC  
.....  
PACKAGE QUANTITY  
50/Case

Adjustable Dial  
Elongated Mask

STANDARD CONNECTOR

REORDER NO. 33072

ADULT  
.....  
PACKAGE QUANTITY  
50/Case

Adjustable Dial  
Elongated Mask

STANDARD CONNECTOR

REORDER NO. 33073

PEDIATRIC  
.....  
PACKAGE QUANTITY  
50/Case



# CO<sub>2</sub> MASKS & ACCESSORIES



Equal  
Length  
CO<sub>2</sub> Sampling  
Line & O<sub>2</sub>  
Tubing

33080 & 33081

**dynarex**  
*Resp-C<sub>2</sub>*

Dynarex Resp-O<sub>2</sub>™ Oxygen CO<sub>2</sub> Masks & Accessories are designed to help monitor patients' exhaled carbon dioxide levels while maintaining the continuous administration of oxygen. The masks feature an elongated shape for added patient comfort and swivel design for optimal flexibility and range of motion.

- Designed to monitor CO<sub>2</sub> levels from the nose and mouth
- Soft, elongated mask provides a secure fit
- Transparent materials allow for clear visualization
- Includes 7' of crush-resistant, star lumen oxygen supply tubing
- 33080 & 33081 feature CO<sub>2</sub> sampling line and O<sub>2</sub> tubing of equal length
- Not made with natural rubber latex or DEHP

### Oxygen CO<sub>2</sub> Elongated Mask

CO<sub>2</sub> SAMPLING LINE WITH FILTER  
7' (2.1 M) OXYGEN TUBE  
MALE/MALE CONNECTOR

**REORDER NO. 33080**

ADULT

PACKAGE QUANTITY  
25/Case

### Oxygen CO<sub>2</sub> Elongated Mask

CO<sub>2</sub> SAMPLING LINE WITH FILTER  
7' (2.1 M) OXYGEN TUBE  
MALE/MALE CONNECTOR

**REORDER NO. 33081**

PEDIATRIC

PACKAGE QUANTITY  
25/Case

Female Connector

Male Connector

FILTER

### Oxygen CO<sub>2</sub> Pressure Line Filter

**REORDER NO. 33087**

ALL AGES

PACKAGE QUANTITY  
50/Case

### Oxygen CO<sub>2</sub> 7' (2.1 m) Sampling Line

FEMALE/FEMALE CONNECTOR

**REORDER NO. 33082**

ALL AGES

PACKAGE QUANTITY  
25/Case

### Oxygen CO<sub>2</sub> 7' (2.1 m) Sampling Line

MALE/MALE CONNECTOR

**REORDER NO. 33083**

ALL AGES

PACKAGE QUANTITY  
25/Case

### Oxygen CO<sub>2</sub> 7' (2.1 m) Sampling Line

MALE/FEMALE CONNECTOR

**REORDER NO. 33084**

ALL AGES

PACKAGE QUANTITY  
25/Case

### Oxygen CO<sub>2</sub> 7' (2.1 m) Sampling Line

MALE/MALE CONNECTOR  
WITH FILTER

**REORDER NO. 33085**

ALL AGES

PACKAGE QUANTITY  
25/Case

### Oxygen CO<sub>2</sub> 7' (2.1 m) Sampling Line

MALE/FEMALE CONNECTOR  
WITH FILTER

**REORDER NO. 33086**


ALL AGES

PACKAGE QUANTITY  
25/Case





# NASAL OXYGEN CANNULAS

## Standard Tip



**dynarex**  
**Resp-O<sub>2</sub>**

Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas with Standard Tips are designed to efficiently deliver supplemental oxygen to respiratory patients. These translucent cannulas are available with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.



#### 4' (1.2 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33100**

ADULT

PACKAGE QUANTITY  
50/Case

#### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33101**

ADULT

PACKAGE QUANTITY  
50/Case

#### 7' (2.1 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33102**

ADULT

PACKAGE QUANTITY  
50/Case

#### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33103**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

#### 7' (2.1 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33104**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

- Non-flared nasal prongs
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP

#### 14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33105**

ADULT

PACKAGE QUANTITY  
50/Case

#### 14' (4.3 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33106**

ADULT

PACKAGE QUANTITY  
50/Case

#### 21' (6.4 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33107**

ADULT

PACKAGE QUANTITY  
25/Case

#### 50' (15.24 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33108**

ADULT

PACKAGE QUANTITY  
15/Case



# NASAL OXYGEN CANNULAS

Flared Tip



dynarex  
**Resp-C<sub>2</sub>**

Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas with Flared Tips are designed to efficiently deliver supplemental oxygen to respiratory patients. The cannulas are available with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Flared-tip nasal prongs help ensure a secure fit
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP

Star Lumen



Standard Connector



7' (2.1 m) Crush  
Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33110**

ADULT

.....  
PACKAGE QUANTITY  
50/Case

7' (2.1 m) Crush  
Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33111**

ADULT

.....  
PACKAGE QUANTITY  
50/Case

14' (4.3 m) Crush  
Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33112**

ADULT

.....  
PACKAGE QUANTITY  
50/Case

21' (6.4 m) Crush  
Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33113**

ADULT

.....  
PACKAGE QUANTITY  
25/Case


Universal Connector





# NASAL OXYGEN CANNULAS

## Cushion Tip



Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas with Cushion Tips are designed to efficiently deliver supplemental oxygen to respiratory patients. The cannulas come with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Lightweight cushioned nasal prongs provide added comfort
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP



### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33120**

NEONATAL

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33121**

INFANT

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33122**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### 4' (1.2 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33123**

ADULT

PACKAGE QUANTITY  
50/Case

### 4' (1.2 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33124**

ADULT

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33125**

ADULT

PACKAGE QUANTITY  
50/Case

### 14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33126**

ADULT

PACKAGE QUANTITY  
50/Case

### 21' (6.4 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33127**

ADULT

PACKAGE QUANTITY  
25/Case

### 25' (7.6 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33128**

ADULT

PACKAGE QUANTITY  
25/Case

### 50' (15.2 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33129**

ADULT

PACKAGE QUANTITY  
15/Case



# NASAL OXYGEN CANNULAS

## Super Soft Touch, Cushion Tip

Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas, Super Soft Touch with Cushion Tips are designed to efficiently deliver supplemental oxygen to respiratory patients. The cannulas come with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Soft vinyl materials and cushioned nasal prongs increase comfort and help minimize the risk of irritation
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP

**dynarex**  
**Resp-O<sub>2</sub>**



### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33140**

NEONATAL

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33141**

INFANT

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33142**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### 4' (1.2 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33143**

ADULT

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33144**

ADULT

PACKAGE QUANTITY  
50/Case

### 14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33145**

ADULT

PACKAGE QUANTITY  
50/Case

### 25' (7.6 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33146**

ADULT

PACKAGE QUANTITY  
25/Case

### 25' (7.6 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33147**

ADULT

PACKAGE QUANTITY  
25/Case



# NASAL OXYGEN CANNULAS

## Cushion Tip with Ultra Soft Ear Foam

Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas with Cushion Tips and Ultra Soft Ear Foam are designed to efficiently deliver supplemental oxygen to respiratory patients. The cannulas come with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Adjustable foam ear cushions and lightweight cushioned nasal prongs increase comfort and help minimize the risk of irritation
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP



7' (2.1 m) Crush  
Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33150**

PEDIATRIC

.....  
PACKAGE QUANTITY  
50/Case

36

7' (2.1 m) Crush  
Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33151**

PEDIATRIC

.....  
PACKAGE QUANTITY  
50/Case

4' (1.2 m) Crush  
Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33152**

ADULT

.....  
PACKAGE QUANTITY  
50/Case

7' (2.1 m) Crush  
Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33153**

ADULT

.....  
PACKAGE QUANTITY  
50/Case

7' (2.1 m) Crush  
Resistant Tubing

UNIVERSAL CONNECTOR


**REORDER NO. 33154**

ADULT

.....  
PACKAGE QUANTITY  
50/Case



Standard Connector



Universal Connector

Star Lumen

## 14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

REORDER NO. 33155

ADULT

PACKAGE QUANTITY  
50/Case

## 25' (7.6 m) Crush Resistant Tubing

STANDARD CONNECTOR

REORDER NO. 33156

ADULT

PACKAGE QUANTITY  
25/Case




dynarex  
Resp-O<sub>2</sub>



# NASAL OXYGEN CANNULAS

Cushion Tip with  
Advantage Ear Foam



### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33170**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33171**

PEDIATRIC

PACKAGE QUANTITY  
50/Case

Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas with Cushion Tips and Advantage Ear Foam are designed to efficiently deliver supplemental oxygen to respiratory patients. These translucent cannulas come with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Adjustable sponge-foam ear wraps and cushioned nasal prongs provide added comfort and help minimize the risk of irritation
- Includes crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP

**33173**  
comes  
with gray  
ear foam

### 4' (1.2 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33172**

ADULT

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33173**

ADULT

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33174**

ADULT

PACKAGE QUANTITY  
50/Case

### 14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33175**

ADULT

PACKAGE QUANTITY  
50/Case

### 25' (7.6 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33176**

ADULT


PACKAGE QUANTITY  
25/Case



# NASAL OXYGEN HIGH FLOW CANNULAS

Cushion Tip





Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen High Flow Cannulas are designed to deliver supplemental oxygen at high flow rates to patients in respiratory distress. The cannulas come with Standard connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Delivers up to 15 liters of oxygen per minute
- Soft, cushioned nasal prongs provide added comfort
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP

### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33160**

ADULT

PACKAGE QUANTITY

25/Case

### 14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33161**

ADULT

PACKAGE QUANTITY

25/Case

COMING  
SOON

### 25' (7.6 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33162**

ADULT

PACKAGE QUANTITY

10/Case



# EtCO<sub>2</sub> NASAL SAMPLING CANNULAS

Equal  
Length EtCO<sub>2</sub>  
Sampling  
Line &  
O<sub>2</sub> Tubing

Dynarex Resp-O<sub>2</sub>™ EtCO<sub>2</sub> Nasal Sampling Cannulas with Cushion Nasal Tips are designed to help monitor exhaled carbon dioxide during respiratory treatment. Each soft, cushioned cannula features a permanent barrier between the nasal prongs so that EtCO<sub>2</sub> can be sampled from one nostril while oxygen continues to flow to the other nostril.

- Permanent barrier between nasal prongs allows for EtCO<sub>2</sub> sampling with uninterrupted oxygen delivery

- Constructed from soft materials for added patient comfort
- Crush-resistant tubing helps prevent occlusion and kinking
- Features EtCO<sub>2</sub> sampling line and O<sub>2</sub> tubing of equal length
- Transparent design allows for clear visualization
- Not made with natural rubber latex or DEHP



Female  
Luer  
Connector



Male  
Luer  
Connector



Standard  
Tubing  
Connector



Universal  
Tubing  
Connector



7' (2.1 m) Tubing

MALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR

REORDER NO. 33180

INFANT

PACKAGE QUANTITY  
25/Case

7' (2.1 m) Tubing

MALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR

REORDER NO. 33181

PEDIATRIC

PACKAGE QUANTITY  
25/Case

7' (2.1 m) Tubing

MALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR

REORDER NO. 33182

ADULT

PACKAGE QUANTITY  
25/Case

**dynarex**  
**Resp-CO<sub>2</sub>**

## 7' (2.1 m) Tubing

FEMALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33183**INFANT  
.....PACKAGE QUANTITY  
25/Case

## 7' (2.1 m) Tubing

FEMALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33184**PEDIATRIC  
.....PACKAGE QUANTITY  
25/Case

## 7' (2.1 m) Tubing

FEMALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33185**ADULT  
.....PACKAGE QUANTITY  
25/Case

## 10' (3.0 m) Tubing

MALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33186**PEDIATRIC  
.....PACKAGE QUANTITY  
25/Case

## 10' (3.0 m) Tubing

MALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33187**ADULT  
.....PACKAGE QUANTITY  
25/Case

## 14' (4.2 m) Tubing

MALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33188**PEDIATRIC  
.....PACKAGE QUANTITY  
25/Case

## 14' (4.2 m) Tubing

MALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33189**ADULT  
.....PACKAGE QUANTITY  
25/Case

## 10' (3.0 m) Tubing

FEMALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33190**PEDIATRIC  
.....PACKAGE QUANTITY  
25/Case

## 10' (3.0 m) Tubing

FEMALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33191**ADULT  
.....PACKAGE QUANTITY  
25/Case

## 14' (4.2 m) Tubing

FEMALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33192**PEDIATRIC  
.....PACKAGE QUANTITY  
25/Case

## 14' (4.2 m) Tubing

FEMALE LUER CONNECTOR  
STANDARD O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33193**ADULT  
.....PACKAGE QUANTITY  
25/Case

## 7' (2.1 m) Tubing

MALE LUER CONNECTOR  
UNIVERSAL O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33194**PEDIATRIC  
.....PACKAGE QUANTITY  
25/Case

## 7' (2.1 m) Tubing

MALE LUER CONNECTOR  
UNIVERSAL O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33195**ADULT  
.....PACKAGE QUANTITY  
25/Case

## 7' (2.1 m) Tubing


FEMALE LUER CONNECTOR  
UNIVERSAL O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33196**ADULT  
.....PACKAGE QUANTITY  
25/Case

## 10' (3.0 m) Tubing

FEMALE LUER CONNECTOR  
UNIVERSAL O<sub>2</sub> TUBING CONNECTOR**REORDER NO. 33197**ADULT  
.....PACKAGE QUANTITY  
25/Case




# OXYGEN SUPPLY TUBING



Dynarex Resp-O<sub>2</sub>™ Oxygen Supply Tubing is designed to extend a patient's reach when accessing a concentrated oxygen source for respiratory therapy. The tubing features a light green color, designed to help improve tube visibility. It comes with a standard connector on one end and either a standard or universal connector on the other.

- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP



### 7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33200**

PACKAGE QUANTITY  
50/Case

### 7' (2.1 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33201**

PACKAGE QUANTITY  
50/Case

### 14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33202**

PACKAGE QUANTITY  
50/Case

### 14' (4.3 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

**REORDER NO. 33203**

PACKAGE QUANTITY  
50/Case

### 21' (6.4 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33204**

PACKAGE QUANTITY  
25/Case

### 25' (7.6 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33205**

PACKAGE QUANTITY  
20/Case

### 30' (9.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33206**

PACKAGE QUANTITY  
20/Case

### 40' (12.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33207**

PACKAGE QUANTITY  
15/Case

### 50' (15.2 m) Crush Resistant Tubing

STANDARD CONNECTOR

**REORDER NO. 33208**

PACKAGE QUANTITY  
15/Case



# OXYGEN CONCENTRATOR & ACCESSORIES

## Oxygen Concentrator

The Dynarex Resp-O<sub>2</sub>™ Oxygen Concentrator is designed to provide supplemental oxygen to help relieve breathing difficulties in respiratory patients. The 5-liter Oxygen Concentrator includes audible alarms, oxygen purity indicator (OPI), and user-friendly flow-meter knob to provide prescribed liters-per-minute. The 37-pound concentrator is easily transportable with accessible top handle and four rolling casters.

- Features a flow range of 0.5-5 LPM
- Audible and visual alarms include Low Flow-Rate, Low/High Pressure, and Power Failure
- Manual reset button helps ensure uninterrupted oxygen delivery
- Designed for quiet operation with low noise output
- Easy-to-read control panel for optimal visibility



## Oxygen Concentrator

5 LITER

REORDER NO. 33950

PACKAGE QUANTITY  
1/Case

**dynarex**  
**Resp-O<sub>2</sub>**

## Sterile Water for Inhalation

Dynarex Resp-O<sub>2</sub>™ Sterile Water supplies are disposable, single-use respiratory therapy bottles filled with water designed for patient inhalation. The water is sterilized to ensure patient safety and to help avoid cross-contamination.

- Reservoir bottles feature built-in oxygen tubing connectors
- Sterile, transparent adaptor body provides optimal visualization
- Audible alarm alerts clinician to air-flow constrictions



**Sterile Water for Inhalation, USP**

**PREFILLED BUBBLE HUMIDIFIER**

**REORDER NO. 34300**

**SIZE**  
350 mL

**PACKAGE QUANTITY**  
20/Case (with adaptor)



**Sterile Water for Inhalation, USP**

**PREFILLED BUBBLE HUMIDIFIER**

**REORDER NO. 34301**

**SIZE**  
550 mL

**PACKAGE QUANTITY**  
12/Case (with adaptor)

## Bubble Humidifier

The Dynarex Resp-O<sub>2</sub>™ Disposable Bubble Humidifier is designed to humidify supplemental oxygen for patients receiving oxygen therapy at home. Ideal for use with all types of oxygen concentrators with either face masks or cannulas, the Bubble Humidifier is intended for oxygen low flow rates up to 6 PSI. The transparent humidifier features clearly defined minimum and maximum lines for easy monitoring.

- Features pressure-relief valve with audible alarm for optimal safety
- Includes PVC diffuser to help minimize noise
- Disposable, single-patient use for optimal hygiene
- Not made with natural rubber latex or DEHP



**Dry Bubble Humidifier**

**400 ML • 6 PSI (HOME USE)**

**REORDER NO. 34221**

**PACKAGE QUANTITY**  
50/Case

## Oxygen Humidifier Connectors

Dynarex Resp-O<sub>2</sub>™ Oxygen Humidifier Connectors are designed to attach dry humidifier bottles or pre-filled bubble humidifiers to continuous-flow oxygen concentrators. They come with either sturdy oxygen tubing on one end and a male threaded grip elbow on the other end, or with star lumen oxygen tubing on one end and molded threaded straight connector on the other end. Connectors are universal and fit most humidifier bottles with a DISS connection.

- Transparent oxygen tubing allows for easy visualization
- Universal fit connects with most humidifier bottles and concentrators
- Ribbed connector provides a secure fit without damaging connection points
- Designed to attach to barb-style concentrator outlets
- Not made with natural rubber latex or DEHP



**Oxygen Humidifier Elbow Connectors with 12" Tubing**

**REORDER NO. 34505**

**PACKAGE QUANTITY**  
10/Case



**Oxygen Humidifier Straight Connectors with 18" Tubing**

**REORDER NO. 34506**

**PACKAGE QUANTITY**  
50/Case



# OXYGEN CYLINDERS





Dynarex Resp-O<sub>2</sub>™ Oxygen Cylinders are premium-grade containers designed to safely and effectively store liquid oxygen for use in medical facilities and homecare settings. The Cylinders are constructed from durable, yet lightweight materials (25% to 40% lighter than steel), making them easy to carry and transport. The brushed aluminum tanks feature a clear coat finish for long-lasting protection. Resp-O<sub>2</sub> Oxygen Cylinders are available in M-6 (B), M-9 (C), M-15 (D) and M-24 (E) models to meet the specific needs of a variety of patients.

- Premium-quality cylinders designed to contain liquid oxygen safely and reliably
- Lightweight construction (25%-40% lighter than steel) allows for easy portability
- Brushed aluminum with clear coat finish for optimal protection
- Standard post-valve oxygen cylinders accept CGA 870-style regulators
- Compliant with 5-year hydrostatic testing
- Ideal for medical settings and homecare use
- Ships empty; prescription for oxygen required for filling



### Oxygen Sleeve Bag

**FOR USE WITH 36391  
OXYGEN D-CYLINDER**

**REORDER NO. 36610**

**PACKAGE QUANTITY**  
1/Case



### M24 Oxygen E-Cylinder

**4.6 L • WITH POST VALVE**

**REORDER NO. 36390**

**PACKAGE QUANTITY**  
1/Case



### M15 Oxygen D-Cylinder

**2.9 L • WITH CGA870 POST VALVE**

**REORDER NO. 36391**

**PACKAGE QUANTITY**  
1/Case



### M9 Oxygen C-Cylinder

**1.7 L • WITH CGA870 POST VALVE**

**REORDER NO. 36392**

**PACKAGE QUANTITY**  
1/Case



### M6 Oxygen B-Cylinder

**1 L • WITH CGA870 POST VALVE**

**REORDER NO. 36393**

**PACKAGE QUANTITY**  
1/Case



# OXYGEN REGULATORS & CYLINDER WRENCH

## Oxygen Regulators

Dynarex Resp-O<sub>2</sub>™ Oxygen Flow Regulators are designed to safely and efficiently administer and regulate the flow of oxygen to respiratory patients at adjustable flow-rate levels. The Oxygen Regulator features an anodized aluminum body with a brass core and brass-encased component parts for optimal safety and reliability. Some select models feature all-brass construction.

- Single-selector knob allows for effortless adjustment
- Features easy-to-read flow-rate dial
- Brass core construction provides increased safety and reliability
- 36350-36361: Comply to CGA 870 standard
- 36370-36374: Comply to CGA 540 standard
- Not made with natural rubber latex or DEHP

## Cylinder Wrench

The Dynarex Resp-O<sub>2</sub>™ Oxygen Cylinder Wrench is designed for use with CGA 870-style cylinders to open and close oxygen-tank valves efficiently and safely. The wrench is constructed from strong, durable plastic and ergonomically crafted to effectively manage a variety of cylinders with minimal exertion.

- Sturdy, break-resistant plastic provides reliable, long-lasting performance.
- Ideal for use in a variety of medical settings, long-term care, and homecare
- Sleek, streamlined design allows for quick and easy adjustments.
- Not made with natural rubber latex or DEHP





36350-36353



36354-36357



36358-36359



36360



36361



36370-36374



Cylinder Wrench

REORDER NO. 36600

PACKAGE QUANTITY  
25/Case

	Reorder No.	Product Features	Flow Range	Dial settings												Qty.
				1	2	3	4	5	6	7	8	9	10	11	12	
CGA 870	36350	BARB OUTLET • BRASS-SLEEVED	0-4 LPM	0	1/32	1/16	1/8	1/4	1/2	3/4	1	1 1/2	2	3	4	5/CASE
	36351	BARB OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/8	1/2	1	1 1/2	2	2 1/2	3	4	5	6	8	5/CASE
	36352	BARB OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE
	36353	BARB OUTLET • BRASS-SLEEVED	0-25 LPM	0	1/2	1	2	3	4	5	6	8	10	15	25	5/CASE
	36354	MINI • BARB OUTLET • BRASS-SLEEVED	0-4 LPM	0	1/32	1/16	1/8	1/4	1/2	3/4	1	1 1/2	2	3	4	5/CASE
	36355	MINI • BARB OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/8	1/2	1	1 1/2	2	2 1/2	3	4	5	6	8	5/CASE
	36356	MINI • BARB OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE
	36357	MINI • BARB OUTLET • BRASS-SLEEVED	0-25 LPM	0	1/2	1	2	3	4	5	6	8	10	15	25	5/CASE
	36358	EXTENDED • DISS OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/8	1/2	1	1 1/2	2	2 1/2	3	4	5	6	8	5/CASE
	36359	EXTENDED • DISS OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE
	36360	BARB OUTLET • ALL BRASS	0-25 LPM	0	1/2	1	2	4	6	8	10	12	15	20	25	5/CASE
	36361	BARB OUTLET • 2 DISS OUTLETS • ALL BRASS	0-25 LPM	0	1/2	1	2	4	6	8	10	12	15	20	25	5/CASE
CGA 540	36370	DISS OUTLET • BRASS-SLEEVED	0-4 LPM	0	1/32	1/16	1/8	1/4	1/2	3/4	1	1 1/2	2	3	4	5/CASE
	36371	DISS OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/2	1	1 1/2	2	2 1/2	3	4	5	6	7	8	5/CASE
	36372	DISS OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE
	36373	MINI • DISS OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/2	1	1 1/2	2	2 1/2	3	4	5	6	7	8	5/CASE
	36374	MINI • DISS OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE



# OXYGEN TANK CARTS AND HOLDERS





Oxygen Tank Cart

**REORDER NO. 36620****PACKAGE QUANTITY**  
1/Case

Dual Oxygen Tank Cart

**REORDER NO. 36621****PACKAGE QUANTITY**  
1/CaseOxygen Tank Holder  
Wheelchair Mount**REORDER NO. 36622****PACKAGE QUANTITY**  
1/Case

Dynarex Resp-O<sub>2</sub>™ Oxygen Tank Carts are designed to safely and reliably carry and support D- and E-sized oxygen cylinders. The sturdy frame is constructed from durable carbon steel with an 80-lb capacity. The cylinder cart is available in single-tank and dual-tank varieties and features a height-adjustable handle for easy customization.

- Durable yet lightweight cart provides optimal portability
- 36620 & 36621: Silver-vein finish is easy to maintain
- 36622: Chrome plating helps to resist corrosion
- Large, composite wheels allow for smoother maneuverability
- 80 lb. weight capacity
- Not made with natural rubber latex or DEHP

**dynarex**  
*Resp-O<sub>2</sub>*™



# OXYGEN THERAPY ACCESSORIES

Dynarex Resp-O<sub>2</sub>™ Oxygen Therapy Accessories are designed to improve the efficacy, convenience, and overall comfort of respiratory equipment. They are not made with natural rubber latex or DEHP.

**Nasal Oxygen Cannula Ear Foam Covers**

- Lightweight foam construction
- Suitable for adult and pediatric patients

**Oxygen Tubing Adaptors**

- Designed to extend adaptability of oxygen tubing

**Inline Water Traps**

- Traps and collects excess moisture from oxygen tubing
- Keeps liquid away from patients' nasal passages
- Tapered fittings ensure a secure connection to the tubing
- Transparent body allows for clear monitoring of fluid

**Portable Oxygen Liter Meter**

- Designed to measure gas-flow from oxygen source
- Convenient pen-style design
- Not to be used as a calibration device



Nasal Oxygen Cannula  
Deluxe Ear Foam Covers

**REORDER NO. 33300**

**PACKAGE QUANTITY**  
50 pairs/Case



Nasal Oxygen Cannula  
Advantage Ear Foam  
Covers

**REORDER NO. 33301**

**PACKAGE QUANTITY**  
50 pairs/Case



Oxygen Tubing  
Inline Water Trap

**REORDER NO. 33302**

**PACKAGE QUANTITY**  
10/Case



Oxygen Tubing 22 mm  
I.D. x 5-7 mm Tapered  
Nipple Adaptor

**REORDER NO. 33303**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing  
5-7 mm End Adaptor

**REORDER NO. 33304**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing  
Swivel Male/Male  
Adaptor

**REORDER NO. 33305**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing  
Y Adaptor

**REORDER NO. 33306**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing  
Nut Nipple Adaptor

**REORDER NO. 33307**

**PACKAGE QUANTITY**  
50/Case



Oxygen Tubing  
Nut Nipple  
Swivel Adaptor

**REORDER NO. 33308**

**PACKAGE QUANTITY**  
50/Case




Portable Oxygen  
Liter Meter

**REORDER NO. 33309**

**PACKAGE QUANTITY**  
15/Case



# CONNECTORS AND ADAPTORS





### Multi-Standard Adaptors

15 MM ID × 22 MM OD  
CONNECTIONS

**REORDER NO. 34510**

**PACKAGE QUANTITY**  
50/Case



### Cuff Connectors

STANDARD 22 MM ID × 22 MM ID  
CONNECTIONS

**REORDER NO. 34511**

**PACKAGE QUANTITY**  
50/Case



### Corrugated Cuff Adaptors

TAPERED TO ACCEPT 22 MM ID  
CUFF TUBING

**REORDER NO. 34512**

**PACKAGE QUANTITY**  
50/Case

Dynarex Resp-O<sub>2</sub>™ Connectors and Adaptors are engineered to provide improved adaptability and connectivity for use in a broad range of respiratory therapy applications—including oxygen therapy, aerosol therapy, and ventilation. These high-quality devices are manufactured from top-grade plastic materials for optimal durability and dependability. Individual models include features such as swivel, cuffed, or barbed designs for added adaptability for multiple connections.

- Designed to optimize connectivity for a variety of respiratory therapy applications
- Engineered from premium plastic materials for optimal reliability
- Translucent design for optimal visibility
- Not made with natural rubber latex or DEHP



### Pressure Line Adaptors

22 MM OD × 22 MM ID CONNECTIONS  
BARBED ELBOW ACCEPTS 5-7 MM ID  
TUBING • ELBOW SWIVELS 360°

**REORDER NO. 34513**

**PACKAGE QUANTITY**  
50/Case



### Straight Adaptors

22 MM OD × 22 MM ID  
WITH OXYGEN STEM

**REORDER NO. 34514**

**PACKAGE QUANTITY**  
50/Case



### Aerosol Tee Connectors

STANDARD 22 MM OD × 15 MM ID ×  
22 MM OD CONNECTIONS

**REORDER NO. 34515**

**PACKAGE QUANTITY**  
50/Case



### Nebulizer Tee Connectors

STANDARD 22 MM OD × 18 MM ID ×  
22 MM OD CONNECTIONS

**REORDER NO. 34516**

**PACKAGE QUANTITY**  
50/Case



# PULSE OXIMETERS



Dynarex Resp-O<sub>2</sub>™ Fingertip Pulse Oximeters are premium-quality, non-invasive devices designed to effectively measure pulse rate (PR) and blood oxygenation levels (SpO<sub>2</sub>). The Pulse Oximeters are designed to fit fingertips comfortably and securely and feature easy-to-read displays for optimal visibility. These compact Oximeters are lightweight and feature a convenient lanyard for easy portability.

- Designed to provide effective measurement of pulse rate and blood oxygenation levels
- Large font display allows for accurate reading
- Constructed from durable materials for reliable performance
- Includes convenient lanyard for easy portability
- 36400: Suitable for all ages
- 36401: Suitable for ages 12 and up; Features Tuya App Bluetooth capabilities
- 36402: Pediatric-sized Oximeter features child-friendly teddy-bear design
- Not made with natural rubber latex or DEHP

36400



36401



Tuya App  
Bluetooth  
capabilities

36402



Child-friendly  
Teddy Bear Design

### Fingertip Pulse Oximeter

ALL AGES

REORDER NO. 36400

PACKAGE QUANTITY  
10/Case

### Fingertip Pulse Oximeter

AGES 12+ • BLUETOOTH CAPABLE

REORDER NO. 36401

PACKAGE QUANTITY  
10/Case

### Fingertip Pulse Oximeter

PEDIATRIC

REORDER NO. 36402

PACKAGE QUANTITY  
10/Case





# STERILE SALINE AND WATER

Dynarex Resp-O<sub>2</sub>™ Sterile Saline and Water supplies are disposable, single-use respiratory therapy bottles/vials filled with water or saline solutions designed for patient inhalation or device irrigation. The water and saline solutions are sterilized to ensure patient safety and to help avoid cross-contamination.

## **Sterile Water for Irrigation, USP & Sterile Normal Saline**

- Suitable for suctioning and irrigation
- Not suitable for injection

## **Sodium Chloride Inhalation Solution, USP – 0.9% Unit Dose Vials**

- Provides humidification of respiratory mucosa and helps prevent dryness
- Tamper-proof, color-coded vials allow for effortless identification
- Graduated markings allow clinicians to control dosage by the drop or milliliter (mL)

## **Prefilled Sterile Water for Inhalation, USP – Nebulizers**

- Transparent container allows for quantity view
- Packaged in a convenient bottle with an easy snap-off cap

## **Sterile Water for Inhalation, USP – Prefilled Bubble Humidifier with Adaptor**

- Reservoir bottles feature built-in oxygen tubing connectors
- Sterile, transparent adaptor body provides optimal visualization
- Audible alarm alerts clinician to air-flow constrictions

**dynarex**  
**Resp-O<sub>2</sub>**

COMING  
SOON

### Sterile Water for Irrigation, USP

NOT FOR INJECTION

**REORDER NO. 34330**

**SIZE**  
100 mL

**PACKAGE QUANTITY**  
48/Case

COMING  
SOON

### Sterile Normal Saline

0.9% SODIUM CHLORIDE  
FOR IRRIGATION, USP  
NOT FOR INJECTION

**REORDER NO. 34331**

**SIZE**  
100 mL

**PACKAGE QUANTITY**  
48/Case



### Sodium Chloride Inhalation Solution, USP

0.9% UNIT DOSE VIAL

**REORDER NO. 34322**

**SIZE**  
3 mL

**PACKAGE QUANTITY**  
1000/Case



### Sodium Chloride Inhalation Solution, USP

0.9% UNIT DOSE VIAL

**REORDER NO. 34323**

**SIZE**  
5 mL

**PACKAGE QUANTITY**  
1000/Case



### Sodium Chloride Inhalation Solution, USP

0.9% UNIT DOSE VIAL

**REORDER NO. 34324**

**SIZE**  
15 mL

**PACKAGE QUANTITY**  
144/Case



### Prefilled Sterile Water for Inhalation, USP

NEBULIZER

**REORDER NO. 34310**

**SIZE**  
500 mL

**PACKAGE QUANTITY**  
12/Case



### Prefilled Sterile Water for Inhalation, USP

NEBULIZER

**REORDER NO. 34311**

**SIZE**  
750 mL

**PACKAGE QUANTITY**  
12/Case



### Prefilled Sterile Water for Inhalation, USP

NEBULIZER

**REORDER NO. 34312**

**SIZE**  
1000 mL

**PACKAGE QUANTITY**  
12/Case



### Sterile Water for Inhalation, USP

PREFILLED BUBBLE HUMIDIFIER

**REORDER NO. 34300**

**SIZE**  
350 mL

**PACKAGE QUANTITY**  
20/Case (with adaptor)



### Sterile Water for Inhalation, USP

PREFILLED BUBBLE HUMIDIFIER

**REORDER NO. 34301**

**SIZE**  
550 mL

**PACKAGE QUANTITY**  
12/Case (with adaptor)



# AEROSOL MASKS

Dynarex Resp-O<sub>2</sub>™ Aerosol Masks are designed for use in aerosol therapy to help administer inhalable treatments and aerosol medications to respiratory patients.

- Transparent vinyl mask features smooth edges
- Mask features an adjustable nose clip and elastic strap
- Can be used with most handheld nebulizers
- Connector accepts 22 mm aerosol corrugated tubing and nebulizer cup
- Not made with natural rubber latex or DEHP

**dynarex**  
**Resp-O<sub>2</sub>**™

STANDARD



ELONGATED



FACE TENT MASK

## Aerosol Face Tent Mask

REORDER NO. 34005

ADULT

PACKAGE QUANTITY  
50/Case

## Aerosol Standard Mask

REORDER NO. 34000

ADULT

PACKAGE QUANTITY  
50/Case

## Aerosol Standard Mask

REORDER NO. 34001

PEDIATRIC

PACKAGE QUANTITY  
50/Case

## Aerosol Standard Mask

REORDER NO. 34002

INFANT

PACKAGE QUANTITY  
50/Case

## Aerosol Elongated Mask

REORDER NO. 34003

ADULT

PACKAGE QUANTITY  
50/Case

## Aerosol Elongated Mask

REORDER NO. 34004


PEDIATRIC

PACKAGE QUANTITY  
50/Case



# TRACHEOSTOMY MASK

with Swivel Tubing Connector



Dynarex Resp-O<sub>2</sub>™ Tracheostomy Masks are designed to help administer aerosol medications and oxygen to tracheostomy patients during respiratory therapy. The collar-style masks feature swivel connectors that allow for access from either side of the patient.

- Collar-style mask features smooth edges
- Non-constricting, adjustable neckband features a snap-closure for easy removal
- Swivel connector accepts 22 mm aerosol corrugated tubing
- Not made with natural rubber latex or DEHP

(5.4 mm) Adaptor  
or Accepts 22 mm  
Aerosol Tubing

**REORDER NO. 34020**

**ADULT**

.....  
**PACKAGE QUANTITY**  
50/Case

One Side Snap • Accepts  
22 mm Aerosol Tubing

**REORDER NO. 34021**

**ADULT**

.....  
**PACKAGE QUANTITY**  
50/Case

One Side Snap • Accepts  
22mm Aerosol Tubing

**REORDER NO. 34022**

**PEDIATRIC**

.....  
**PACKAGE QUANTITY**  
50/Case



**dynarex**  
**Resp-C<sub>2</sub>**



# AEROSOL DRAINAGE BAGS

## 750 cc Bag

## Accepts Standard 22 mm Corrugated Tubing

Dynarex Resp-O<sub>2</sub>™ Aerosol Drainage Bags are designed to prevent accumulated liquid and condensation from reaching patients' airways during vent/aerosol treatment.

- Traps and contains accumulated liquid and condensation
- Self-locking plastic chain secures bag to stable objects
- Not made with natural rubber latex or DEHP



**Y-ADAPTOR**



**T-ADAPTOR**

Y-Adaptor

**REORDER NO. 34200**

**PACKAGE QUANTITY**  
50/Case

T-Adaptor

**REORDER NO. 34201**

**PACKAGE QUANTITY**  
50/Case

# DISPOSABLE WATER TRAPS

## Y-Cap



The Dynarex Resp-O<sub>2</sub>™ Disposable Water Trap is designed to seal off excess moisture and gasses during oxygen therapy and prevent them from backing up into tubing. The Water Trap features a threaded collection jar and self-sealing Y-shaped lid to ensure that circuits remain clear and unobstructed.

- Designed to trap liquids and help prevent obstructions in corrugated/circuit tubing
- Features threaded jar and Y-sealing lid to ensure uninterrupted flow
- Constructed from high-quality plastic materials for optimal reliability
- Not made with natural rubber latex or DEHP

For use with  
22 mm Tubing

**REORDER NO. 34202**

**ADULT**

.....  
**60 CC**  
.....

**PACKAGE QUANTITY**  
50/Case

For use with  
15 mm Tubing

**REORDER NO. 34203**

**PEDIATRIC**

.....  
**60 CC**  
.....

**PACKAGE QUANTITY**  
50/Case


**dynarex**  
**Resp-C<sub>2</sub>**™



# SMALL VOLUME NEBULIZER 6 CC CUPS

Dynarex Resp-O<sub>2</sub>™ Small Volume Nebulizer Cups are an integral part of the nebulizing process, in which liquid medication is transformed into an inhaled mist to aid in respiratory therapy treatments. The cups feature an anti-spill design that allows patients to receive treatments upright or at an angle.

- Effectively delivers medication at angles up to 90°
- Features a threaded cap to help prevent leaks
- Disposable cup capacity: 6 mL
- Suitable for hand-held and table-top nebulizers
- Not made with natural rubber latex or DEHP



**REORDER NO. 34100**

**ACCESSORIES**  
6 cc cup only

.....  
**ALL AGES**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34101**

**ACCESSORIES**  
7' (2.1 m) Tubing  
Standard Connector

.....  
**ALL AGES**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34102**

**ACCESSORIES**  
7' (2.1 m) Tubing  
Standard Connector  
Spring Loaded T  
(22 mm OD × 22 mm OD)

.....  
**ALL AGES**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34103**

**ACCESSORIES**  
7' (2.1 m) Tubing  
Standard Connector  
Spring Loaded T  
(22 mm ID × 22 mm OD)

.....  
**ALL AGES**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34104**

**ACCESSORIES**  
7' (2.1 m) Tubing  
Standard Connector  
T-Piece  
Mouthpiece

.....  
**ALL AGES**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34105**

**ACCESSORIES**  
7' (2.1 m) Tubing  
Standard Connector  
T-Piece  
Mouthpiece  
6" Flex Tube

.....  
**ALL AGES**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34106**

**ACCESSORIES**  
7' (2.1 m) Tubing  
Universal Connector  
T-Piece  
Mouthpiece  
6" Flex Tube

.....  
**ALL AGES**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34107**

**ACCESSORIES**  
7' (2.1 m) Tubing  
Standard Connector  
45 Degree Pacifier Elbow

.....  
**INFANT/TODDLER**  
.....

**PACKAGE QUANTITY**  
12/Case

**REORDER NO. 34120**

**ACCESSORIES**  
Aerosol Elongated Mask  
7' (2.1 m) Tubing  
Standard Connector

.....  
**ADULT**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34121**

**ACCESSORIES**  
Aerosol Elongated Mask  
7' (2.1 m) Tubing  
Universal Connector

.....  
**ADULT**  
.....

**PACKAGE QUANTITY**  
50/Case

**REORDER NO. 34122**

**ACCESSORIES**  
Aerosol Elongated Mask  
7' (2.1 m) Tubing  
Standard Connector

.....  
**PEDIATRIC**  
.....

**PACKAGE QUANTITY**  
50/Case


**REORDER NO. 34123**

**ACCESSORIES**  
Aerosol Elongated Mask  
7' (2.1 m) Tubing  
Universal Connector


.....  
**PEDIATRIC**  
.....

**PACKAGE QUANTITY**  
50/Case

Aerosol  
Elongated Mask







# SPRING LOADED T-ADAPTORS


**34500**
**22 mm OD × 22 mm OD**
**REORDER NO. 34500**
**PACKAGE QUANTITY**  
30/Case

**34501**
**22 mm ID × 22 mm OD**
**REORDER NO. 34501**
**PACKAGE QUANTITY**  
30/Case

**34502**
**15 mm ID × 15 mm OD**
**REORDER NO. 34502**
**PACKAGE QUANTITY**  
30/Case

Dynarex Resp-O<sub>2</sub>™ Spring Loaded T-Adaptors are designed to be used in conjunction with a small volume nebulizer (SVN) to help deliver aerosol medication to respiratory vent patients. The T-Adaptor features a spring-loaded valve that opens and closes the nebulizer circuit during insertion and removal.

- Eliminates the need to break the ventilator circuit between treatments
- Ventless design prevents medication leakage
- Accepts most standard aerosol tubing/ventilator circuit
- Not made with natural rubber latex or DEHP



# PEDIATRIC COMPRESSOR NEBULIZERS



**dynarex**  
*Resp-C<sub>2</sub>* 

Dynarex Resp-O<sub>2</sub>™ Pediatric Compressor Nebulizers are lightweight medical devices, specifically designed for children, that allow for the easy inhalation of prescribed medication. These adorable, child-friendly designs help caretakers manage children's respiratory disorders easily and effectively.

- Small Volume Nebulizer (SVN)
- Converts liquid medication into an aerosol mist for easy inhalation
- Quiet operation levels
- Includes reusable nebulizer cup
- Not made with natural rubber latex or DEHP



**ACCESSORIES**

<b>INPUT</b>	AC 120 V, 60 HZ
<b>EXTREME PRESSURE</b>	>205 kpa
<b>POWER CONSUMPTION</b>	≤200 VA
<b>MEDICATION CAPACITY</b>	10 mL
<b>PARTICLE SIZE</b>	≤5 μm & 60%
<b>MMAD</b>	2.0 μm
<b>AVERAGE NEBULIZATION RATE</b>	≥0.2 mL/min
<b>NOISE LEVEL</b>	≤ 60 dB
<b>COMPRESSOR PRESSURE RANGE</b>	30-58 psi
<b>OPERATION PRESSURE RANGE</b>	7-10 psi
<b>LITER FLOW RANGE</b>	6-10 lpm
<b>ACCESSORIES</b>	Infant Mask, Pediatric Mask, Reusable Nebulizer Cup, Mouthpiece, T-Piece, 7' (2.1 m) Tube, 5 Filters



**Carry Bag**

**REORDER NO. 34406**

**DIMENSIONS**  
8.3" × 6.5" × 5.5"

**PACKAGE QUANTITY**  
4/Case



**My Bear**

**REORDER NO. 34402**

**DIMENSIONS**  
8.3" × 6.5" × 5.5"

**PACKAGE QUANTITY**  
4/Case



**My Doggy**

**REORDER NO. 34403**

**DIMENSIONS**  
7.1" × 6.3" × 9.1"

**PACKAGE QUANTITY**  
4/Case



**My Piggy**

**REORDER NO. 34404**

**DIMENSIONS**  
8.3" × 6.5" × 5.5"

**PACKAGE QUANTITY**  
4/Case



**My Penguin**

**REORDER NO. 34405**

**DIMENSIONS**  
6.7" × 5.1" × 6.5"

**PACKAGE QUANTITY**  
4/Case



# COMPRESSOR NEBULIZERS

Dynarex Resp-O<sub>2</sub>™ Compressor Nebulizers feature compact, ergonomic designs that allow for the easy inhalation of prescribed medication. These nebulizers help users to efficiently manage respiratory ailments.

- Small Volume Nebulizer (SVN)
- Converts prescribed medication into an aerosol mist for easy inhalation
- Quiet operation levels
- Includes reusable nebulizer cup
- Not made with natural rubber latex or DEHP

<b>INPUT</b>	AC 120 V, 60 HZ
<b>EXTREME PRESSURE</b>	>205 kpa
<b>POWER CONSUMPTION</b>	≤200 VA
<b>MEDICATION CAPACITY</b>	10 mL
<b>PARTICLE SIZE</b>	≤5 μm & 60%
<b>MMAD</b>	2.0 μm
<b>AVERAGE NEBULIZATION RATE</b>	≥0.2 mL/min
<b>NOISE LEVEL</b>	≤ 60 dB
<b>COMPRESSOR PRESSURE RANGE</b>	30-58 psi
<b>OPERATION PRESSURE RANGE</b>	14.5-27 psi (1 bar-1.9 bar)
<b>LITER FLOW RANGE</b>	16-23.5 lpm
<b>ACCESSORIES</b>	Adult Mask, Pediatric Mask, Reusable Nebulizer Cup, Mouthpiece, T-Piece, 7' (2.1 m) Tube, 5 Filters



**Elite**

**REORDER NO. 34400**

**DIMENSIONS**  
29 × 18 × 11 cm

**PACKAGE QUANTITY**  
4/Case



**Compact**

**REORDER NO. 34401**

**DIMENSIONS**  
15.5 × 13 × 8.3 cm

**PACKAGE QUANTITY**  
6/Case



**Compact (International)**

**REORDER NO. 34408**

**DIMENSIONS**  
15.5 × 13 × 8.3 cm

**PACKAGE QUANTITY**  
6/Case



**dynarex**  
**Resp-C<sub>2</sub>**



# 50 PSI COMPRESSOR & ACCESSORIES

## 50 PSI Compressor

The Dynarex Resp-O<sub>2</sub>™ 50 PSI Compressor is a top-quality respiratory device, engineered to efficiently deliver high-pressure compressed air on a continuous basis. The powerful Compressor is designed to administer humidified air or aerosolized medication into patients' airways to aid in the treatment of respiratory conditions. The Compressor is easy to clean and features a sleek, lightweight design for convenient portability. The Compressor features an oil-less piston motor for quiet, low-vibration operation. The Dynarex Resp-O<sub>2</sub>™ Compressor includes a U-shaped DISS adaptor for optimal connectivity.

- Powerful compressor delivers high-pressure compressed air for extended periods of time
- Designed to efficiently deliver humidified air or aerosolized medication
- Sleek, lightweight design (14.9 lb.) allows for optimal portability
- Oil-less piston motor provides quiet, low-vibration operation (noise level of <49.5 dB)
- U-shaped DISS adaptor for optimal connectivity
- Includes two extra air intake filters for added convenience
- Easy to wipe clean for optimal hygiene
- Ideal for use with heated nebulizers, humidifier/nebulizers, and mist tents



50 PSI Compressor

**REORDER NO. 33952**

**PACKAGE QUANTITY**  
1/Case

**dynarex**  
**Resp-O<sub>2</sub>**

## Large Volume Nebulizer

The Dynarex Resp-O<sub>2</sub>™ Large Volume Nebulizer is engineered for safe and effective administration of humidified gas in concentrations ranging from 28 percent to 98 percent. The Nebulizer bottle is constructed from high-quality, break-resistant plastic for optimal durability. The empty Nebulizer Bottle is designed for use with clinician-provided water, offering an economical alternative to pre-filled humidifier bottles. The Nebulizer Bottle is equipped with a Venturi-style entrainment ring with an incremental dial for setting oxygen concentrations. The transparent bottle features clear marking for easy monitoring of water levels.

- Designed to efficiently deliver humidified gas for respiratory therapy
- Constructed from break-resistant plastic for optimal reliability
- Venturi-style entrainment ring with incremental dial for adjusting oxygen concentrations
- FIO<sub>2</sub> settings are adjustable from 28% to 98%
- Transparent plastic design provides optimal visibility of contents
- Clearly marked water levels for quick and easy monitoring
- Not made with natural rubber latex or DEHP



Large Volume Nebulizer

**REORDER NO. 34225**

**PACKAGE QUANTITY**  
24/Case

## Corrugated Tubing

Dynarex Resp-O<sub>2</sub>™ Corrugated Tubing is designed to help administer nebulized medication during aerosol therapy and intermittent positive pressure breathing (IPPB) therapy. This high-quality respiratory tubing is designed from premium vinyl materials for optimal durability. The flexible, lightweight tubing features tapered cuffs on both ends for secure attachment to standard 22mm connections. The disposable Tubing is conveniently segmented every six inches to allow for accurate cutting and sizing. The Tubing features a blue tint design to allow for easy monitoring of condensation build-up. The Tubing is available in 6-foot and 100-foot lengths to suit a variety of needs.

- Designed to help administer aerosol and IPPB applications
- Crafted from high-quality vinyl materials for optimal durability
- Flexible and lightweight for convenient handling
- Tapered cuffs attach to standard 22mm adapters for a secure fit
- Segmented every 6 inches for accurate cutting and sizing
- Blue tint allows for easy visualization of condensation build-up
- Disposable, single-patient use for optimal hygiene
- Available in 6-foot and 100-foot lengths to suit a variety of needs



100 ft. (30.5 m)  
Corrugated Tubing

**WITH STANDARD 22MM ADAPTORS**  
**BLUE SEGMENTED ROLL**

**REORDER NO. 34211**

**PACKAGE QUANTITY**  
1/Case



6ft (1.8 m)  
Corrugated Tubing

**WITH STANDARD 22MM ADAPTORS**  
**BLUE SEGMENTED ROLL**

**REORDER NO. 34213**

**PACKAGE QUANTITY**  
50/Case



# SINGLE SUCTION CATHETERS

Graduated with Control Valve  
Coiled Packed

Dynarex Resp-O<sub>2</sub>™ Single Suction Catheters are designed to safely and effectively remove secretions from the respiratory tract of tracheostomy patients. The soft, flexible Catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. The Catheter includes a vented cap port for enhanced suction control and two staggered eyelet openings to increase suctioning effectiveness and reduce patient trauma.

- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows for greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Coil packed for compact storage
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP

**6 Fr**
**Pediatric**
**REORDER NO. 32000**
**PACKAGE QUANTITY**  
50/Case

**8 Fr**
**Pediatric**
**REORDER NO. 32001**
**PACKAGE QUANTITY**  
50/Case

**10 Fr**
**Adult/Pediatric**
**REORDER NO. 32002**
**PACKAGE QUANTITY**  
50/Case

**12 Fr**
**Adult**
**REORDER NO. 32003**
**PACKAGE QUANTITY**  
50/Case

**14 Fr**
**Adult**
**REORDER NO. 32004**
**PACKAGE QUANTITY**  
50/Case

**16 Fr**
**Adult**
**REORDER NO. 32005**
**PACKAGE QUANTITY**  
50/Case

**18 Fr**
**Adult**
**REORDER NO. 32006**
**PACKAGE QUANTITY**  
50/Case



# SUCTION CATHETER KITS

## CONTENTS:

Graduated Catheter with Control Valve  
Pair of Powder-Free Vinyl Exam Gloves  
Pop-Up Basin

Dynarex Resp-O<sub>2</sub>™ Single Suction Catheter Kits are designed to safely and effectively remove secretions from the respiratory tract of tracheostomy patients. The soft, flexible catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. The catheter features a vented cap control valve for enhanced suction control and two staggered eyelet openings to increase suctioning effectiveness and reduce patient trauma. Each Kit includes a pop-up basin to hold irrigation fluids and a pair of powder-free vinyl exam gloves.

- Kit includes pop-up basin and exam gloves for added convenience
- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows for greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP

6 Fr

Pediatric

REORDER NO. 32050

PACKAGE QUANTITY  
50/Case

8 Fr

Pediatric

REORDER NO. 32051

PACKAGE QUANTITY  
50/Case

10 Fr

Adult/Pediatric

REORDER NO. 32052

PACKAGE QUANTITY  
50/Case

12 Fr

Adult

REORDER NO. 32053

PACKAGE QUANTITY  
50/Case

14 Fr

Adult

REORDER NO. 32054

PACKAGE QUANTITY  
50/Case

16 Fr

Adult

REORDER NO. 32055

PACKAGE QUANTITY  
50/Case

18 Fr

Adult

REORDER NO. 32056

PACKAGE QUANTITY  
50/Case



# SUCTION CATHETER MINI-KITS

## CONTENTS:

Graduated Catheter with Control Valve  
Pair of Powder-Free Vinyl Exam Gloves  
Mini-Tray



Dynarex Resp-O<sub>2</sub>™ Suction Catheter Mini-Kits are designed to safely and effectively remove secretions from the respiratory tract of tracheostomy patients. The soft, flexible catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. The catheter includes a vented cap control valve for enhanced suction control and features two eyelet openings for increased suctioning effectiveness with minimal patient trauma. Each Kit also features a mini-tray to hold irrigation solutions and a pair of powder-free vinyl exam gloves.

- Convenient Mini-Kit includes a pair of vinyl exam gloves and mini-tray to hold irrigation fluids
- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP


**8 Fr**
**Pediatric**
**REORDER NO. 32080**
**PACKAGE QUANTITY**  
100/Case

**10 Fr**
**Adult/Pediatric**
**REORDER NO. 32081**
**PACKAGE QUANTITY**  
100/Case

**12 Fr**
**Adult**
**REORDER NO. 32082**
**PACKAGE QUANTITY**  
100/Case

**14 Fr**
**Adult**
**REORDER NO. 32083**
**PACKAGE QUANTITY**  
100/Case

**16 Fr**
**Adult**
**REORDER NO. 32084**
**PACKAGE QUANTITY**  
100/Case



# SLEEVED SUCTION CATHETER KITS

## CONTENTS:

Graduated Sleeved Catheter with Control Valve  
Pop-Up Basin

Dynarex Resp-O<sub>2</sub>™ Sleeved Suction Catheter Kits are designed to safely and efficiently remove secretions from the respiratory tracts of tracheostomy patients. The soft, flexible catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. Each catheter is contained in a plastic sleeve to provide a protective layer between the patient and clinician. Each Kit includes a pop-up basin to hold irrigation fluids.



- Protective sleeve helps prevent the risk of cross-contamination
- Convenient Kit includes a pop-up basin to hold irrigation fluids
- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows for greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth number markings
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP


**8 Fr**
**Pediatric**
**REORDER NO. 32090**
**PACKAGE QUANTITY**  
50/Case

**10 Fr**
**Adult/Pediatric**
**REORDER NO. 32091**
**PACKAGE QUANTITY**  
50/Case

**12 Fr**
**Adult**
**REORDER NO. 32092**
**PACKAGE QUANTITY**  
50/Case

**14 Fr**
**Adult**
**REORDER NO. 32093**
**PACKAGE QUANTITY**  
50/Case



# SUCTION CONNECTING TUBING

## Non-Conductive • with Male Connector

Dynarex Resp-O<sub>2</sub>™ Suction Connecting Tubing is designed to provide secure connections to both portable and in-wall suction sources. This non-conductive Tubing features female connectors on both ends and includes an optional male (straw) adapter for additional attachments.

- Designed to connect suction sources to suction catheters and yankauers
- Constructed from transparent material for clear visualization during suctioning
- Durable, crush-resistant tubing
- Single-use and disposable
- Not made with natural rubber latex or DEHP



Combo Pack:  $\frac{3}{16}$ "  $\times$  6'  
&  $\frac{3}{16}$ "  $\times$  18" Tubes

**REORDER NO. 32100**

ALL AGES

NON-STERILE

PACKAGE QUANTITY  
50/Case

$\frac{3}{16}$ "  $\times$  6' Tube

**REORDER NO. 32101**

ALL AGES

NON-STERILE

PACKAGE QUANTITY  
50/Case

Combo Pack:  $\frac{3}{16}$ "  $\times$  6'  
&  $\frac{3}{16}$ "  $\times$  18" Tubes

**REORDER NO. 32102**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

$\frac{3}{16}$ "  $\times$  18" Tube

**REORDER NO. 32103**

ALL AGES

STERILE

PACKAGE QUANTITY  
100/Case

$\frac{3}{16}$ "  $\times$  6' Tube

**REORDER NO. 32104**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

$\frac{3}{16}$ "  $\times$  10' Tube

**REORDER NO. 32105**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

$\frac{3}{16}$ "  $\times$  12' Tube

**REORDER NO. 32106**

ALL AGES

STERILE

PACKAGE QUANTITY  
20/Case

$\frac{3}{16}$ "  $\times$  20' Tube

**REORDER NO. 32107**

ALL AGES

STERILE

PACKAGE QUANTITY  
20/Case

$\frac{1}{4}$ "  $\times$  6' Tube

**REORDER NO. 32108**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

$\frac{1}{4}$ "  $\times$  10' Tube

**REORDER NO. 32109**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

$\frac{1}{4}$ "  $\times$  12' Tube

**REORDER NO. 32110**

ALL AGES

STERILE

PACKAGE QUANTITY  
20/Case

$\frac{1}{4}$ "  $\times$  20' Tube

**REORDER NO. 32111**

ALL AGES

STERILE

PACKAGE QUANTITY  
20/Case

$\frac{9}{32}$ "  $\times$  6' Tube

**REORDER NO. 32112**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case



# YANKAUER SUCTION HANDLES





Dynarex Resp-O<sub>2</sub>™ Yankauer Suction Handles are designed to safely and effectively remove particulates and fluids from patients' airways. The curved, one-piece Handle is constructed from transparent materials for clear visualization, and features a slip-resistant, ribbed design to ensure a secure connection to suction tubing.

- Designed to clear airways to help improve ventilation and oxygenation
- Constructed from shatter-resistant plastic materials
- Available in vented, non-vented, bulb-tip, and straight-tip designs
- Sterile, single-use, and disposable for optimal hygiene
- Not made with natural rubber latex or DEHP



#### Bulb Tip Vented

**REORDER NO. 32120**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

#### Bulb Tip Non-Vented

**REORDER NO. 32121**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

#### Straight Tip Vented

**REORDER NO. 32122**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

#### Straight Tip Non-Vented

**REORDER NO. 32123**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

#### Bulb Tip • Non-Vented with Pre-Attached 3/16" × 6' Suction Tubing

**REORDER NO. 32125**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

#### Bulb Tip • Non-Vented with Pre-Attached 1/4" × 6' Suction Tubing

**REORDER NO. 32126**

ALL AGES

STERILE

PACKAGE QUANTITY  
50/Case

#### Bulb Tip • Non-Vented with Pre-Attached 3/16" × 10' Suction Tubing

**REORDER NO. 32127**

ALL AGES

STERILE

PACKAGE QUANTITY  
20/Case

#### Bulb Tip • Non-Vented with Pre-Attached 1/4" × 10' Suction Tubing

**REORDER NO. 32128**

ALL AGES

STERILE

PACKAGE QUANTITY  
20/Case



# DISPOSABLE SUCTION CANISTERS

Dynarex Resp-O<sub>2</sub>™ Disposable Suction Canisters are designed to safely and effectively collect and contain aspirated secretions from patients in medical settings. The Canisters feature an ultrasonically welded shut-off valve to help prevent overflow and cross-contamination. These versatile Canisters are adaptable for use with floor, cabinet, or wall-mounted units. Each disposable, single-use Canister comes with a convenient lid equipped with a pour spout, vacuum, and patient port.

- Constructed from high-quality, durable polymer materials
- Port caps are engineered to create a snug seal to prevent vacuum loss
- Transparent canister features easy-to-read gradations for accurate monitoring
- Not made with natural rubber latex or DEHP



## Disposable Suction Canister 800 cc

FLOAT VALVE SHUTOFF

**REORDER NO. 32200**

**PACKAGE QUANTITY**  
24/Case

## Disposable Suction Canister 800 cc

FLOAT VALVE SHUTOFF  
3/16" x 6' AND 3/16" x 18" TUBING  
WITH MALE CONNECTOR

**REORDER NO. 32201**

**PACKAGE QUANTITY**  
12/Case

## Disposable Suction Canister 1200 cc

FLOAT VALVE SHUTOFF

**REORDER NO. 32220**

**PACKAGE QUANTITY**  
40/Case

## Disposable Suction Canister 1200 cc

FLOAT VALVE SHUTOFF  
3/16" x 6' AND 3/16" x 18" TUBING  
WITH MALE CONNECTOR

**REORDER NO. 32221**

**PACKAGE QUANTITY**  
30/Case

## Disposable EMS Suction Canister 1200 cc

FLOAT VALVE SHUTOFF

**REORDER NO. 32224**

**PACKAGE QUANTITY**  
40/Case

# OUTER SUCTION CANISTERS & LINERS

## Disposable Suction Liners

Dynarex Resp-O<sub>2</sub>™ Disposable Suction Liners are designed to effectively and safely collect and contain aspirated fluids and materials. Each high-quality, durable Liner features a patient-port elbow adaptor that allows quick and easy hookup and helps prevent tubing from kinking. The Liner is ultrasonically sealed to the lid—so no assembly is needed. Suction Liners require the use of a reusable outer canister.

- Ultrasonically welded shut-off valve helps prevent overflow and contamination
- Lid comes sealed and pre-assembled with pour spout, vacuum, and patient ports
- All port caps are engineered to create an air-tight snuggle seal
- Not made with natural rubber latex or DEHP



## Reusable Outer Suction Canisters

Dynarex Resp-O<sub>2</sub>™ Reusable Outer Suction Canisters are designed for holding inner-liners to safely collect and contain secretions and bodily fluids during suctioning. The rigid reusable Canister is precision-crafted from high-quality polymer materials. The Canister features a built-in bracket engineered for easy wall-mounting, and pre-attached vacuum tubing for secure connection with inner-liners.

- Features built-in adjustable on-off valve
- Clearly marked critical measure numbering
- Designed for use with disposable inner-liners
- Not made with natural rubber latex or DEHP



Reusable Outer  
Suction Canisters  
1000 cc

WALL MOUNT

REORDER NO. 32250

PACKAGE QUANTITY  
10/Case

Reusable Outer  
Suction Canisters  
1500 cc

WALL MOUNT

REORDER NO. 32251

PACKAGE QUANTITY  
10/Case

Reusable Outer  
Suction Canisters  
3000 cc

WALL MOUNT

REORDER NO. 32252

PACKAGE QUANTITY  
10/Case

Disposable  
Suction Liner  
1000 cc

REORDER NO. 32260

PACKAGE QUANTITY  
50/Case

Disposable  
Suction Liner  
1500 cc

REORDER NO. 32261

PACKAGE QUANTITY  
50/Case

Disposable  
Suction Liner  
3000 cc

REORDER NO. 32262

PACKAGE QUANTITY  
50/Case



# HOMECARE SUCTION UNIT

## Portable • Battery Operated

The Resp-O<sub>2</sub>™ Homecare Portable Suction Unit is a premium-quality suction unit designed to provide convenient, portable suctioning for use at home or during travel. The Unit comes with a rechargeable battery as well as AC/DC power cord for optimal versatility. The Suction Unit features a maximum vacuum of 600 ±20 mm Hg and a flow range of 22 ±10% LPM. The Suction Unit comes fully equipped with suction bacteria filter, elbow connector, elbow tube, and 6' suction tubing. The Resp-O<sub>2</sub> Homecare Suction Unit is equipped with an easy-to-read gauge to provide precise measurements of vacuum levels when needed for pediatric patients.

- Designed for convenient, effective suctioning for home use and travel
- Features compact, lightweight design (3.6 lb.) for optimal portability
- Maximum vacuum of 600 ±20 mm Hg and a flow range of 22 ±10% LPM
- Equipped with suction bacteria filter, elbow connector, elbow tube, and 6' suction tubing
- Includes an 800 cc disposable suction canister
- Easy-to-read gauge for precise measurement
- Rechargeable Lithium Ion battery provides up to 50 minutes of continuous operation
- Includes AC/DC power cord for added versatility



Homecare Suction Unit

**REORDER NO. 32305**

**PACKAGE QUANTITY**  
1/Case

**dynarex**  
**Resp-O<sub>2</sub>**

# HOME CARE 800 cc SUCTION DISPOSABLE CANISTER KIT

The Dynarex Resp-O<sub>2</sub>™ Homecare Suction Disposable Canister Kit is designed to collect and contain aspirated secretions from patients in medical settings safely and effectively. Each Canister Kit comes equipped with the necessary components for suction therapy. The canister features an ultrasonically welded shut-off valve to help prevent overflow and cross-contamination. These versatile Kits are adaptable for use with most portable suction units. The Kit includes an 800 cc canister, 6' suction tubing, cuttable 16" silicone tubing, elbow connector, and a suction bacteria filter.



- Convenient kit includes necessary accessories for suction therapy
- Constructed from high-quality, durable polymer materials
- Port caps are engineered to create a snug seal to prevent vacuum loss
- Ultrasonically welded shut-off valve helps prevent overflow and cross-contamination
- Transparent canister features easy-to-read gradations for accurate monitoring
- Kit comes equipped with 800 cc canister, 6' suction tubing, cuttable 16" silicone tubing, elbow connector, and suction bacteria filter
- Compatible with most portable suction units
- Not made with natural rubber latex or DEHP

## Suction Disposable Canister

**800 CC  
FLOAT VALVE FULL KIT**

**REORDER NO. 32202**

**PACKAGE QUANTITY**  
1 Kit/Box • 8 Kits/Case



# SUCTION ACCESSORIES



## Sodium Chloride Inhalation Solution, USP – 0.9% Unit Dose Vials

Dynarex Resp-O<sub>2</sub>™ Sodium Chloride Inhalation Solution, USP, 0.9% Unit Dose Vials are disposable, single-use respiratory therapy vials filled with a sterile, normal saline solution (0.9% sodium chloride). The vials are designed to be used with a variety of tracheal suction systems for inhalation, nebulization, and irrigation. The saline solution is sterilized to ensure patient safety and to help avoid cross-contamination.

- Designed for safe and efficient irrigation during suction procedures
- Twist-and-pull top allows for easy opening
- Features graduations in 0.5 mL increments for accurate measurement
- Tapered tip facilitates use with a variety of tracheal suction systems for optimal adaptability



### Sodium Chloride Inhalation Solution, USP

0.9% UNIT DOSE VIAL • 3 ML

**REORDER NO. 34322**

**PACKAGE QUANTITY**  
1000/Case



### Sodium Chloride Inhalation Solution, USP

0.9% UNIT DOSE VIAL • 5 ML

**REORDER NO. 34323**

**PACKAGE QUANTITY**  
1000/Case



### Sodium Chloride Inhalation Solution, USP

0.9% UNIT DOSE VIAL • 15 ML

**REORDER NO. 34324**

**PACKAGE QUANTITY**  
144/Case

**dynarex**  
**Resp-O<sub>2</sub>**

## Fluid Solidifiers

Dynarex Resp-O<sub>2</sub>™ Fluid Solidifiers are designed to be used in conjunction with suction canisters to provide safe and effective containment of biohazardous liquids. The Fluid Solidifier is formulated to quickly turn secretions into a gel, greatly reducing the risk of exposure to contaminated bodily fluids and aerosolized bloodborne pathogens. The Fluid Solidifier allows for the safe transport of infectious waste and eliminates the risk of spills and sprays.

- Formulated to safely encapsulate biohazardous liquids
- Screw-cap bottle allows for easy dispensing and secure storage
- Helps reduce turnover time between procedures by quickly solidifying waste

## Mucus Specimen Trap

The Dynarex Resp-O<sub>2</sub>™ Mucus Specimen Trap is designed for use with suction devices to safely collect specimens and trach cultures for testing and analysis. The transparent 40 cc container features clearly marked gradations in 1 cc increments for easy visualization of contents. The sterile trap includes a disposable adapter to help prevent cross-contamination, and a leak-resistant screw-on cap to help maintain a sterile pouring surface.

- Features heavy-walled, collapse-resistant tubing
- Includes threaded cap for secure seal
- Packed in a convenient, easy-peel envelope
- Includes patient label for easy identification
- Sterile trap is disposable and designed for single use
- Not made with natural rubber latex or DEHP

## Suction Canister Bracket and Suction Canister Ring

Dynarex Resp-O<sub>2</sub>™ Suction Canister Brackets and Rings are designed for use with disposable and reusable canisters for easy attachment and wall mounting. Brackets are constructed from top-grade plastic materials and Rings are made from high-quality stainless steel.



### Fluid Solidifier

**ABSORBS UP TO 500 CC  
OF LIQUID MEDICAL WASTE**

**REORDER NO. 32405**

**PACKAGE QUANTITY**  
96/Case



### Fluid Solidifier

**ABSORBS UP TO 1500 CC  
OF LIQUID MEDICAL WASTE**

**REORDER NO. 32406**

**PACKAGE QUANTITY**  
96/Case



### Fluid Solidifier

**ABSORBS UP TO 3000 CC  
OF LIQUID MEDICAL WASTE**

**REORDER NO. 32407**

**PACKAGE QUANTITY**  
96/Case



### Mucus Specimen Trap

**LEAK-PROOF TRANSPORT CAP  
40 CC**

**REORDER NO. 32270**

**PACKAGE QUANTITY**  
50/Case



### Suction Canister Bracket

**REORDER NO. 32408**

**PACKAGE QUANTITY**  
10/Case



### Suction Canister Ring

**REORDER NO. 32409**

**PACKAGE QUANTITY**  
10/Case



# SUCTION UNITS

Vacuum Levels up to 560 mmHg





32300

Suction Unit with  
Suction Cup Legs

**REORDER NO. 32300**

**PACKAGE QUANTITY**  
1/Case



32301

Suction Unit  
with Base

**REORDER NO. 32301**


**PACKAGE QUANTITY**  
1/Case

Dynarex Resp-O<sub>2</sub>™ Suction Units are premium-grade medical devices designed to safely and effectively remove excess secretions and fluids from the airways of respiratory patients. These powerful suction units are easily portable and feature high-quality vacuum regulators and anti-vibration vacuum gauge settings to ensure quiet and efficient operation.

- Designed to provide vacuum levels up to 560 mmHg
- Fiber filter guards against bacteria and contaminants
- Clearly marked gauge allows for easy pressure monitoring
- Oil-less lubrication provides consistent functionality
- Includes an 800 cc disposable suction canister, tubing and bacterial filter



# TRACHEOSTOMY CARE KITS



**dynarex**  
**Resp-C<sub>2</sub>**

Dynarex Resp-O<sub>2</sub>™ Tracheostomy Care Kits are designed to provide the essential components needed for the preparation, cleaning, and care of tracheostomy tubes. These time-saving Kits are disposable, portable, and compact for easy use and storage. The components in these fully sealed, single-use Kits are sterile and come packed in sequence of use for maximum efficiency.

- Fully sealed and sterile package and contents
- Components are packed in order of procedural use
- Disposable, single-patient use
- Not made with natural rubber latex or DEHP

REORDER NUMBER	35000	35001	35002	35003	35004	35005	35006	35008	35009
TRAY TYPE	ONE COMPARTMENT	THREE COMPARTMENT	TWO COMPARTMENT	TWO COMPARTMENT	TWO COMPARTMENT	THREE COMPARTMENT SALINE	THREE COMPARTMENT SALINE	ONE COMPARTMENT	MINI
REMOVABLE PLASTIC BASIN				1	1				
PAIR OF VINYL STRETCH GLOVES, L	1	1	1	1			1		
PAIR OF NITRILE GLOVES, M					1	1			
3-PLY WATERPROOF DRAPE, 45 × 34 CM	1	1	1	1	1	1	1	1	
NON-WOVEN GAUZE, 4" × 4"	4	4	4	4	4	4	4	4	2
NON-WOVEN DRESSING, 4" × 4"	1	1	1	1	1	1	1	1	2
COTTON TIP APPLICATOR, 5"	2	2	2	2	2	2	2	2	2
TWILL TAPE, 34"	1	1	1	1	1	1	1	1	
PIPE CLEANER, 5.5"	2	2	2	2	2	2	2	2	2
TRACHEOSTOMY BRUSH	1	1	1	1	1	1	1	1	
14 FR. SUCTION CATHETER			1						
STERILE SALINE 0.9%, 100 ML						1			
HYDROGEN PEROXIDE 1 OZ. PACK						2	2		1
PACKAGING QUANTITY	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	20/CASE	40/CASE



# TRACHEOSTOMY CARE KIT COMPONENTS



One-Compartment Tray



Mini-Tray



Removable Plastic Basin



Two-Compartment Tray



Three-Compartment Tray



Three-Compartment Saline Tray



Pair of Vinyl Stretch Gloves, Size L



Pair of Nitrile Gloves, Size M



14 Fr. Suction Catheter



Twill Tape, 34"



3-Ply Waterproof Drape, 45 × 34 cm



Non-Woven Gauze, 4" × 4"



Hydrogen Peroxide 1 oz. Pack



Sterile Saline 0.9%, 100 mL




Cotton Tip Applicator, 5"



Tracheostomy Brush



Pipe Cleaner, 5.5"



Tracheostomy Dressing,  
Non-Woven, 4" × 4"



# TRACHEOSTOMY SPONGES





Dynarex Resp-O<sub>2</sub>™ Tracheostomy Drain Sponges are premium-quality sterile dressings designed with a convenient pre-cut T-shaped notch for a comfortable fit around drains, catheters, chest tubes, IVs, and trach tubes. These reliable, highly absorbent pads are designed to help keep patients dry and comfortable.

- Pre-cut notch saves time and eliminates need for scissors
- Packed in peel-down pouches for sanitary handling
- Pouches come in open trays for easy access
- Sterile, not made with natural rubber latex

### Tracheostomy Drain Sponges

NON-WOVEN • 4" × 4" • 6-PLY

**REORDER NO. 35050**

#### PACKAGE QUANTITY

2 Sponges/Pouch  
25 Pouches/Tray  
12 Trays/Case

### Tracheostomy Gauze Sponges

100% COTTON • 4" × 4" • 12-PLY

**REORDER NO. 35051**

#### PACKAGE QUANTITY

2 Sponges/Pouch  
50 Pouches/Tray  
6 Trays/Case



# COMFORT-FOAM TRACH TUBE HOLDERS

Dynarex Resp-O<sub>2</sub>™ Comfort-Foam Trach Tube Holders are designed to stabilize and hold tracheostomy tubes securely in place to ensure uninterrupted airflow to trach patients. The Trach Tube Holder fastens securely with hook-and-loop closures and features soft, Comfort-Foam padding to help keep skin dry and promote optimal patient comfort.



## Anti-Disconnect • 1-Piece

SMALL

**REORDER NO. 35024**

NEONATAL &amp; INFANT

NECK

7"-9"/18-23 cm

PACKAGE QUANTITY

12/Box • 5 Boxes/Case

## Anti-Disconnect • 1-Piece

MEDIUM

**REORDER NO. 35025**

ADOLESCENT &amp; ADULT

NECK

9"-14"/23-35 cm

PACKAGE QUANTITY

12/Box • 5 Boxes/Case

## Anti-Disconnect • 1-Piece

LARGE

**REORDER NO. 35026**

ADULT

NECK

13"-19"/33-48 cm

PACKAGE QUANTITY

12/Box • 5 Boxes/Case

- Helps to stabilize trach tube and reduce risk of dislodgement
- Comfort-Foam cushioning helps to minimize skin irritation
- Features high-quality, USA-made materials
- Hook-and-loop closures are easy to adjust
- Not made with natural rubber latex or DEHP

## 1-Piece

SMALL

**REORDER NO. 35020**

NEONATAL &amp; INFANT

NECK

7"-9"/18-23 cm

PACKAGE QUANTITY

12/Box • 5 Boxes/Case

## 1-Piece

MEDIUM

**REORDER NO. 35021**

ADOLESCENT &amp; ADULT

NECK

9"-14"/23-35 cm

PACKAGE QUANTITY

12/Box • 5 Boxes/Case

## 1-Piece

LARGE

**REORDER NO. 35022**

ADULT

NECK

13"-19"/33-48 cm

PACKAGE QUANTITY

12/Box • 5 Boxes/Case

## 1-Piece

X-LARGE

**REORDER NO. 35023**

ADULT &amp; BARIATRIC

NECK

16-24"/41-61 cm

PACKAGE QUANTITY

12/Box • 5 Boxes/Case

## 2-Piece

SMALL

**REORDER NO. 35027**

NEONATAL &amp; INFANT

MAXIMUM NECK SIZE

9"/23 cm

PACKAGE QUANTITY

10/Box • 5 Boxes/Case

## 2-Piece

MEDIUM

**REORDER NO. 35028**

PEDIATRIC

MAXIMUM NECK SIZE

18"/46 cm

PACKAGE QUANTITY

10/Box • 5 Boxes/Case

## 2-Piece

LARGE

**REORDER NO. 35029**

ADULT

MAXIMUM NECK SIZE

19.5"/50 cm

PACKAGE QUANTITY

10/Box • 5 Boxes/Case

## 2-Piece

X-LARGE

**REORDER NO. 35030**

BARIATRIC

MAXIMUM NECK SIZE

28"/71 cm

PACKAGE QUANTITY

10/Case



# BREATHABLE-FOAM TRACH TUBE HOLDERS

Dynarex Resp-O<sub>2</sub>™ Breathable-Foam Trach Tube Holders are designed to stabilize and securely hold tracheostomy tubes in place to help ensure uninterrupted airflow to trach patients. The Trach Tube Holder is designed with soft, open-cell foam to aerate skin for optimal patient comfort.

- Helps to stabilize trach tube and reduce risk of dislodgement
- Breathable foam material helps keep skin dry and comfortable
- Features high-quality, USA-made materials
- Adjustable hook-and-loop straps provide secure fit
- Soft-edged outer ties help prevent irritation
- Not made with natural rubber latex or DEHP



## 1-Piece

### SMALL

**REORDER NO. 35031**

#### ADOLESCENT & ADULT

#### NECK

9"-14"/23-35 cm

#### PACKAGE QUANTITY

20/Box • 6 Boxes/Case

## 1-Piece

### MEDIUM

**REORDER NO. 35032**

#### ADULT

#### NECK

11"-17"/28-43 cm

#### PACKAGE QUANTITY

20/Box • 6 Boxes/Case

## 1-Piece

### LARGE

**REORDER NO. 35033**

#### ADULT

#### NECK

13-19"/33-48 cm

#### PACKAGE QUANTITY

20/Box • 6 Boxes/Case

## 2-Piece

### LARGE

**REORDER NO. 35034**

#### ADULT

#### MAXIMUM NECK SIZE

19.5"/50 cm

#### PACKAGE QUANTITY

20/Box • 6 Boxes/Case

# SOFT-FOAM TRACH TUBE HOLDERS



Dynarex Resp-O<sub>2</sub>™ Soft-Foam Trach Tube Holders are designed to stabilize and securely hold tracheostomy tubes in place to help ensure uninterrupted airflow to trach patients. The Trach Tube Holder is constructed from soft, fabric-covered foam that is both strong and lightweight for optimal durability and patient comfort. The Trach Tube Holder features a convenient, one-piece design that is easy to adjust with hook-and-loop straps.

- Helps to stabilize trach tube and reduce risk of dislodgement
- Lightweight, fabric-covered foam provides optimal comfort
- Features high-quality, USA-made materials
- Adjustable hook-and-loop straps provide secure fit
- Not made with natural rubber latex or DEHP

## 1-Piece

### SMALL

**REORDER NO. 35035**

### INFANT

LENGTH  
13.22"

PACKAGE QUANTITY  
100/Case

## 1-Piece

### MEDIUM

**REORDER NO. 35036**

### PEDIATRIC

LENGTH  
15.25"

PACKAGE QUANTITY  
100/Case

## 1-Piece

### LARGE

**REORDER NO. 35037**

### ADULT

LENGTH  
22"

PACKAGE QUANTITY  
100/Case



# MANUAL PULMONARY RESUSCITATORS (MPR)

Dynarex Resp-O<sub>2</sub>™ Manual Pulmonary Resuscitators (MPR) are premium-grade medical devices designed for resuscitation and life support measures performed by qualified medical professionals. The lightweight, single-patient MPR bag features a textured design to provide a non-slip grip and comes pre-assembled with 10 feet of oxygen tubing and soft cushioned facial mask.

- MPR/intake valve ensures one-way air-flow
- Swivel port rotates 360° for optimal freedom of movement
- Translucent components allow for improved visualization
- Disposable, single-patient use
- Not made with natural rubber latex or DEHP



**dynarex**  
**Resp-O<sub>2</sub>**™



Mask • Reservoir Bag  
Oxygen Tubing

**REORDER NO. 36000**

ADULT

.....  
PACKAGE QUANTITY  
6/Case

Mask • Reservoir Tubing  
Oxygen Tubing

**REORDER NO. 36001**

ADULT

.....  
PACKAGE QUANTITY  
6/Case

Mask • Reservoir Bag  
Oxygen Tubing • PEEP Valve

**REORDER NO. 36002**

ADULT

.....  
PACKAGE QUANTITY  
6/Case

Mask • Reservoir Bag  
Oxygen Tubing

**REORDER NO. 36003**

PEDIATRIC

.....  
PACKAGE QUANTITY  
6/Case

Mask • Reservoir Tubing  
Oxygen Tubing

**REORDER NO. 36004**

PEDIATRIC

.....  
PACKAGE QUANTITY  
6/Case

Mask • Reservoir Bag  
Oxygen Tubing

**REORDER NO. 36005**

INFANT

.....  
PACKAGE QUANTITY  
6/Case

Mask • Reservoir Tubing  
Oxygen Tubing

**REORDER NO. 36006**

INFANT

.....  
PACKAGE QUANTITY  
6/Case

Mask • Reservoir Bag  
Oxygen Tubing • PEEP Valve

**REORDER NO. 36007**

PEDIATRIC

.....  
PACKAGE QUANTITY  
6/Case



# AIR CUSHION MASKS WITH VALVE



Dynarex Resp-O<sub>2</sub>™ Air Cushion Masks are premium-grade masks designed for the administration of anesthesia and/or for use in conjunction with manual pulmonary resuscitation (MPR) bags. The Cushion Mask features a facial-contour design and inflatable check valve to provide a secure seal with minimal pressure.

- Features ergonomic shape for a snug, comfortable fit
- Constructed from transparent materials for clear visualization
- Air-filled cushioned face plate provides maximum comfort
- Disposable, single-use
- Not made with natural rubber latex or DEHP



Size #0

No Hook

**REORDER NO. 36030**

PREMIE/NEONATAL

 .....  
**PACKAGE QUANTITY**  
 20/Case


Size #1

No Hook

**REORDER NO. 36031**

NEONATE/INFANT

 .....  
**PACKAGE QUANTITY**  
 20/Case


Size #2

Yellow Hook

**REORDER NO. 36032**

TODDLER

 .....  
**PACKAGE QUANTITY**  
 20/Case


Size #3

Green Hook

**REORDER NO. 36033**

CHILD/SMALL ADULT

 .....  
**PACKAGE QUANTITY**  
 20/Case


Size #4

Red Hook

**REORDER NO. 36034**

MEDIUM ADULT

 .....  
**PACKAGE QUANTITY**  
 20/Case


Size #5

Blue Hook

**REORDER NO. 36035**

LARGE ADULT

 .....  
**PACKAGE QUANTITY**  
 20/Case



# PULMONARY FUNCTION MEASUREMENT

## Spirometers

Dynarex Resp-O<sub>2</sub>™ Incentive Spirometers are hand-held devices designed to help improve lung function and minimize the risk of post-operative respiratory complications. The compact, durable spirometer helps motivate patients to strengthen their breathing muscles, particularly when recovering from prolonged anesthesia, surgery, or extended bed rest. The device features clear, easy-to-read markings to provide visual feedback to help patients monitor their progress.

- Adjustable indicator provides target guidelines to indicate sufficient lung exercise
- Flexible mouthpiece allows for easy adjustment and optimal comfort
- Features dual-sided, clearly marked calibrations for effortless visualization
- Compact design provides enhanced efficiency and ease of handling
- Not made with natural rubber latex or DEHP

## Peak Flow Meters

Dynarex Resp-O<sub>2</sub>™ Peak Flow Meters are hand-held monitoring devices that measure peak expiratory flow (PEF). This lightweight, durable Meter is designed to help patients track daily changes in their breathing patterns. The Dynarex Resp-O<sub>2</sub>™ Peak Flow Meter is available in Adult size (standard range 50-800 L/min) and Pediatric size (low range 50-300 L/min).

- Color-coded indicators help to track daily progress
- Easy-grip handle for steady handling
- Not made with natural rubber latex or DEHP
- Includes peak flow rate diary to track daily progress
- Disposable mouthpieces are available for added protection against bacteria and cross-contamination



### Volumetric Incentive Spirometer

2500 mL

**REORDER NO. 36100**

ALL AGES

PACKAGE QUANTITY  
12/Case



### Volumetric Incentive Spirometer

5000 mL

**REORDER NO. 36101**

ALL AGES

PACKAGE QUANTITY  
12/Case



### Incentive Spirometer

TRI-FLOW CHAMBER

**REORDER NO. 36102**

ALL AGES

PACKAGE QUANTITY  
12/Case



### Peak Flow Meter

LOW RANGE 50-300 L/MIN

**REORDER NO. 36110**

PEDIATRIC

PACKAGE QUANTITY  
10/Case



### Peak Flow Meter

FULL RANGE 50-800 L/MIN

**REORDER NO. 36111**

ADULT

PACKAGE QUANTITY  
10/Case

### PFM Nose Clips

The Dynarex Resp-O<sub>2</sub>™ PFM Nose Clip is designed to help prevent nasal air leakage during spirometry and pulmonary function procedures. The premium-quality Nose Clip features a sturdy spring mechanism for a secure hold and soft foam padding for a comfortable fit.

- Soft foam pads help ensure patient comfort
- Spring mechanism keeps nose clip firmly in place
- Disposable, single-patient use for optimal hygiene
- Not made with natural rubber latex or DEHP



### Durable Peak Flow Mouthpiece

WITH ONE-WAY VALVE

**REORDER NO. 36116**

ADULT

PACKAGE QUANTITY  
100/Case



### Cardboard Peak Flow Mouthpiece

DISPOSABLE

**REORDER NO. 36117**

PEDIATRIC

PACKAGE QUANTITY  
100/Case



### Cardboard Peak Flow Mouthpiece

DISPOSABLE

**REORDER NO. 36118**

ADULT

PACKAGE QUANTITY  
100/Case



### Peak Flow Meter Nose Clip

DISPOSABLE

**REORDER NO. 36119**

PEDIATRIC/ADULT

PACKAGE QUANTITY  
25/Bag, 4 Bags/Case



# ENDOTRACHEAL TUBES

## Uncuffed

Dynarex Resp-O<sub>2</sub>™ Endotracheal Tubes are high-quality plastic tubes designed to safely and effectively deliver higher concentrations of oxygen and medications directly into patients' lungs.

- High-quality, durable plastic tubing helps secure and maintain patients' airways
- Designed to deliver oxygen in higher concentrations directly into patients' lungs
- Can be used to safely administer medications directly into the lungs
- Features smoothly crafted Murphy eye holes to help protect vocal cords
- Available cuffed or uncuffed, and with or without intubation stylette
- Stylettes are engineered to help support and guide the tube through patients' airways
- Sterile, not made with natural rubber latex



Uncuffed • 2.0 mm

**REORDER NO. 36220**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 2.5 mm

**REORDER NO. 36221**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 3.0 mm

**REORDER NO. 36222**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 3.5 mm

**REORDER NO. 36223**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 4.0 mm

**REORDER NO. 36224**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 4.5 mm

**REORDER NO. 36225**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 5.0 mm

**REORDER NO. 36226**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 5.5 mm

**REORDER NO. 36227**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 2.5 mm  
with Stylette

**REORDER NO. 36228**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 3.0 mm  
with Stylette

**REORDER NO. 36229**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 3.5 mm  
with Stylette

**REORDER NO. 36230**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 4.0 mm  
with Stylette

**REORDER NO. 36231**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 4.5 mm  
with Stylette

**REORDER NO. 36232**

STERILE

PACKAGE QUANTITY  
10/Box

Uncuffed • 5.0 mm  
with Stylette

**REORDER NO. 36233**

STERILE

PACKAGE QUANTITY  
10/Box




# ENDOTRACHEAL TUBES

## Cuffed

Dynarex Resp-O<sub>2</sub>™ Endotracheal Tubes are high-quality plastic tubes designed to safely and effectively deliver higher concentrations of oxygen and medications directly into patients' lungs.

- High-quality, durable plastic tubing helps secure and maintain patients' airways
- Designed to deliver oxygen in higher concentrations directly into patients' lungs
- Can be used to safely administer medications directly into the lungs
- Features smoothly crafted Murphy eye holes to help protect vocal cords
- Available cuffed or uncuffed, and with or without intubation stylette
- Stylettes are engineered to help support and guide the tube through patients' airways
- Sterile, not made with natural rubber latex



Cuffed • 4.0 mm

**REORDER NO. 36234**

**STERILE**

**PACKAGE QUANTITY**  
10/Box

Cuffed • 4.5 mm

**REORDER NO. 36235**

**STERILE**

**PACKAGE QUANTITY**  
10/Box

Cuffed • 5.0 mm

**REORDER NO. 36236**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 5.5 mm

**REORDER NO. 36237**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 6.0 mm

**REORDER NO. 36238**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 6.5 mm

**REORDER NO. 36239**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 7.0 mm

**REORDER NO. 36240**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 7.5 mm

**REORDER NO. 36241**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 8.0 mm

**REORDER NO. 36242**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 8.5 mm

**REORDER NO. 36243**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 9.0 mm

**REORDER NO. 36244**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 4.0 mm  
with Stylette**REORDER NO. 36245**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 4.5 mm  
with Stylette**REORDER NO. 36246**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 5.0 mm  
with Stylette**REORDER NO. 36247**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 5.5 mm  
with Stylette**REORDER NO. 36248**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 6.0 mm  
with Stylette**REORDER NO. 36249**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 6.5 mm  
with Stylette**REORDER NO. 36250**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 7.0 mm  
with Stylette**REORDER NO. 36251**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 7.5 mm  
with Stylette**REORDER NO. 36252**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 8.0 mm  
with Stylette**REORDER NO. 36253**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 8.5 mm  
with Stylette**REORDER NO. 36254**

STERILE

PACKAGE QUANTITY

10/Box

Cuffed • 9.0 mm  
with Stylette**REORDER NO. 36255**


STERILE

PACKAGE QUANTITY

10/Box



# ENDOTRACHEAL TUBE INTRODUCER (BOUGIE)



Dynarex Resp-O<sub>2</sub>™ Endotracheal Tube Introducers (Bougies) are premium medical devices designed to effectively guide and support endotracheal tubes during difficult intubation procedures. The sterile Introducer is disposable and comes individually wrapped for optimal hygiene.

- Designed to guide and support endotracheal tubes into patients' airways
- Ideal for difficult intubations and "blind" intubations
- Crafted from high-quality plastic materials for optimal durability
- Features clearly visible depth markings for improved safety
- Available in Child and Adult sizes
- Disposable and individually wrapped for optimal hygiene
- Not made with natural rubber latex

**10 Fr**

**Bougie**

**REORDER NO. 36260**

**STERILE**

**PACKAGE QUANTITY**  
10/Box

**15 Fr**

**Bougie**

**REORDER NO. 36261**

**STERILE**

**PACKAGE QUANTITY**  
10/Box

# INTUBATING STYLETTES



Dynarex Resp-O<sub>2</sub>™ Intubating Stylettes are designed to facilitate safe and effective endotracheal intubations. The Stylettes feature a malleable aluminum wire encased in a plastic sheath to help reduce friction for smoother insertion and removal of endotracheal tubes.

- Designed to facilitate safe and efficient endotracheal intubations
- Malleable aluminum wire helps to guide endotracheal tube
- Plastic sheath helps reduce friction for smoother insertion and removal
- Plastic covering extends beyond aluminum wire to minimize trauma
- Available in four sizes to fit most adult and pediatric tracheal tubes
- Ideal for hospital and EMS use
- Not made with natural rubber latex

**6 Fr**

Intubating Stylettes

**REORDER NO. 36265**

STERILE

PACKAGE QUANTITY  
20/Box

**10 Fr**

Intubating Stylettes

**REORDER NO. 36266**

STERILE

PACKAGE QUANTITY  
20/Box

**12 Fr**

Intubating Stylettes

**REORDER NO. 36267**

STERILE

PACKAGE QUANTITY  
20/Box

**14 Fr**

Intubating Stylettes

**REORDER NO. 36268**

STERILE

PACKAGE QUANTITY  
20/Box



# NASOPHARYNGEAL AIRWAYS & KITS



Dynarex Resp-O<sub>2</sub>™ Nasopharyngeal Airways are designed for nasal insertion to help create and maintain an open and secure airway. The Nasopharyngeal Airway features a fixed flange end for safer insertion and internal striations for optimal suctioning efficiency. The Nasopharyngeal Airway Kits include an assortment of airways in multiple sizes for optimal versatility and packets of lubrication jelly for added convenience.

**5.0 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36275**

20 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

**5.5 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36276**

22 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

**6.0 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36277**

24 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

- Designed to be inserted into nasal passage
- Fixed flange design limits the depth of insertion for safe, secure placement
- Constructed from soft materials for enhanced comfort
- Features rounded tip for gentle insertion
- Strong, flexible tubing resists collapse and kinking
- Not made from natural rubber latex

Kit of 6

**REORDER NO. 36270**

AIRWAY SIZES INCLUDED (MM)  
5.5, 6.0, 6.5, 7.0, 7.5, 8.0  
.....  
6 DYNALUBE PACKETS  
.....  
PACKAGE QUANTITY  
1 Kit/Box

Kit of 9

**REORDER NO. 36271**

AIRWAY SIZES INCLUDED (MM)  
5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0  
.....  
9 DYNALUBE PACKETS  
.....  
PACKAGE QUANTITY  
1 Kit/Box

**6.5 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36278**

26 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

**7.0 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36279**

28 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

**7.5 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36280**

30 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

**8.0 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36281**

32 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

**8.5 mm**

Nasopharyngeal  
Airway

**REORDER NO. 36282**

34 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box

**9.0 mm**

Nasopharyngeal  
Airway


**REORDER NO. 36283**

36 FR  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Box



# LARYNGEAL MASK AIRWAYS

## Silicone



Dynarex Resp-O<sub>2</sub>™ Laryngeal Mask Airways (LMA) are premium-quality respiratory devices designed to help open and maintain airways in unconscious patients. The Laryngeal Mask Airway features a cushioned inflatable cuff with reinforced wire coil to provide a direct, sealed route to the lungs. The Laryngeal Mask Airway is ideal for use by anesthesiologists administering general anesthesia during surgery and other medical procedures. The LMA also provides a simpler alternative to endotracheal intubation in emergency situations, allowing for relatively easy insertion and blind intubation, even by novice operators.

- Inflatable cuff fills oral pharynx to help protect airway
- Soft materials and contoured design help reduce throat irritation
- Ideal for use during the administration of anesthesia
- Allows for relatively easy placement and blind intubation
- Translucent tubing and depth marks allow for accurate placement
- Reinforced with wire coil for optimal support and reduced risk of kinking
- Single-use mask is sterile and x-ray detectable

1.0 mm  
Reinforced

REORDER NO. 36285

≤5 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
5/Case

1.5 mm  
Reinforced

REORDER NO. 36286

5-10 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
5/Case

2.0 mm  
Reinforced

REORDER NO. 36288

10-20 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
5/Case

2.5 mm  
Reinforced

REORDER NO. 36289

20-30 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
5/Case

3.0 mm  
Reinforced

REORDER NO. 36290

30-50 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
5/Case

4.0 mm  
Reinforced

REORDER NO. 36291

50-70 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
5/Case

5.0 mm  
Reinforced

REORDER NO. 36292

70-100 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
5/Case

1.0 mm  
Non-Reinforced

REORDER NO. 36293

5 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Case

1.5 mm  
Non-Reinforced

REORDER NO. 36294

5-10 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Case

2.0 mm  
Non-Reinforced

REORDER NO. 36295

10-20 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Case

2.5 mm  
Non-Reinforced

REORDER NO. 36296

20-30 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Case

3.0 mm  
Non-Reinforced

REORDER NO. 36297

30-50 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Case

4.0 mm  
Non-Reinforced

REORDER NO. 36298

50-70 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Case

5.0 mm  
Non-Reinforced

REORDER NO. 36299

70-100 KG  
.....  
STERILE  
.....  
PACKAGE QUANTITY  
10/Case



# ORAL AIRWAYS (BERMAN TYPE)

Dynarex Resp-O<sub>2</sub>™ Oral Berman Airways are designed to aid the breathing process in unconscious patients by preventing the tongue from blocking the posterior pharyngeal wall. The Berman Airway features a dual-channel design with perforated midway openings to facilitate the insertion of suction catheters without obstructing the air passage. The Berman Airway features a smooth, curved design that contours to the shape of the pharynx for improved comfort. The Airway features a sturdy bite block for added stability and color-coding for easy size identification.

- Constructed from premium-quality polyethylene for optimal reliability
- Slotted sides allow for smooth catheter insertion without obstructing air passageway
- Midway openings help maintain passageway even during airway shifts or partial occlusions
- Smooth, non-brittle finish allows for improved flexibility and durability
- Anatomically curved design helps provide improved comfort and easy insertion
- Color-coded for easy size identification
- Individually packaged for optimal hygiene
- Not made with natural rubber latex or DEHP



Kit of 6

**REORDER NO. 36312**

**SIZES INCLUDED**

40, 60, 80, 90, 100, 110 mm

**PACKAGE QUANTITY**

1/Case

40 mm

REORDER NO. 36300

INFANT

PACKAGE QUANTITY

100/Box

40 mm

REORDER NO. 36301

INFANT

PACKAGE QUANTITY

24/Box

60 mm

REORDER NO. 36302

CHILD

PACKAGE QUANTITY

100/Box

60 mm

REORDER NO. 36303

CHILD

PACKAGE QUANTITY

24/Box

80 mm

REORDER NO. 36304

ADULT SMALL

PACKAGE QUANTITY

100/Box

80 mm

REORDER NO. 36305

ADULT SMALL

PACKAGE QUANTITY

24/Box

90 mm

REORDER NO. 36306

ADULT MEDIUM

PACKAGE QUANTITY

100/Box

90 mm

REORDER NO. 36307

ADULT MEDIUM

PACKAGE QUANTITY

24/Box

100 mm

REORDER NO. 36308

ADULT LARGE

PACKAGE QUANTITY

100/Box

100 mm

REORDER NO. 36309

ADULT LARGE

PACKAGE QUANTITY

24/Box

110 mm

REORDER NO. 36310

ADULT X-LARGE

PACKAGE QUANTITY

100/Box

110 mm

REORDER NO. 36311

ADULT X-LARGE

PACKAGE QUANTITY

24/Box

# DISPOSABLE GUEDEL AIRWAYS



**dynarex**  
**Resp-C<sub>2</sub>**



The Dynarex Resp-O<sub>2</sub>™ Guedel Airway is a premium-quality respiratory device designed to help open and maintain airways in unconscious patients. The Guedel Airway's curved construction helps to ensure an uninterrupted breathing process while minimizing patient discomfort. The Airways feature a durable bite block for added stability and are color-coded for easy size identification.

- Constructed from premium-quality polypropylene for optimal durability
- Designed with smooth, curved edges for improved comfort
- Reinforced bite block provides optimal stability
- Transparent design allows for effortless visualization
- Color coding allows for immediate size identification
- Individually packaged for optimal hygiene
- Not made with natural rubber latex or DEHP


**40 mm**
**REORDER NO. 36315**
**NEWBORN**
**PACKAGE QUANTITY**  
10/Box

**50 mm**
**REORDER NO. 36316**
**INFANT**
**PACKAGE QUANTITY**  
10/Box

**60 mm**
**REORDER NO. 36317**
**TODDLER**
**PACKAGE QUANTITY**  
10/Box

**70 mm**
**REORDER NO. 36318**
**CHILD**
**PACKAGE QUANTITY**  
10/Box

**80 mm**
**REORDER NO. 36319**
**ADULT SMALL**
**PACKAGE QUANTITY**  
10/Box

**90 mm**
**REORDER NO. 36320**
**ADULT MEDIUM**
**PACKAGE QUANTITY**  
10/Box

**100 mm**
**REORDER NO. 36321**
**ADULT LARGE**
**PACKAGE QUANTITY**  
10/Box

**110 mm**
**REORDER NO. 36322**
**ADULT X-LARGE**
**PACKAGE QUANTITY**  
10/Box

**120 mm**
**REORDER NO. 36323**
**ADULT XX-LARGE**
**PACKAGE QUANTITY**  
10/Box

# CPR PRODUCTS

Dynarex Resp-O<sub>2</sub>™ CPR Shields and Masks are designed to provide a physical barrier between the responder and patient, reducing the risk of cross-contamination during CPR. The CPR Shield features a barrier filter and one-way valve that prevents patient blowback to the caregiver throughout administration.

- Serves as a protective barrier to prevent direct contact during CPR
- Features a one-way-valve to keep blowback from reaching the rescuer
- Designed from moldable plastic to effectively seal around the nose and mouth
- Individually poly bagged and constructed from disposable materials for optimal hygiene
- Not made with natural rubber latex or DEHP
- 36021: Features convenient key-ring pouch
- 36022: All-in-one kit includes mask, optional oxygen port, vinyl gloves, and hard case



CPR Face Shield

WITH ONE-WAY VALVE  
AND BARRIER FILTER

**REORDER NO. 36020**

**PACKAGE QUANTITY**  
100/Case



CPR Face Shield in  
Soft Case

WITH ONE-WAY VALVE  
AND BARRIER FILTER

**REORDER NO. 36021**

**PACKAGE QUANTITY**  
100/Case



CPR Rescue Mask Kit

WITH ONE-WAY VALVE, BARRIER FILTER  
AND PAIR OF VINYL EXAM GLOVES

**REORDER NO. 36022**

**PACKAGE QUANTITY**  
25/Case







dynarex  
*Resp-C<sub>2</sub>*™



# ADDITIONAL INFO

**Disclaimer:** Dynarex will make every effort to inform our customers of any changes made to a specific product, catalog copy or marketing collateral materials. Dynarex cannot be held responsible for these changes. Dynarex reserves the right to make changes to the technical specs, dimensions, look, or any other product change at any time without prior notice except as required by the US FDA Code of Federal Regulations for changes made to medical devices, OTC drugs and other foreign governmental regulations.

Go to [www.dynarex.com](http://www.dynarex.com) for the latest product information. Customers should consult with their Customer Care Representative before placing a purchase order. Please note: All products in the Dynarex catalog may not be licensed in all countries. Customers are hereby directed to consult with their Customer Care Representative before placing a purchase order.

Product sizes may vary slightly based on production tolerances using the metric system and are stated for marketing purposes only.

Go to [www.dynarex.com](http://www.dynarex.com) to find which products are labeled with the California Proposition 65 warning to indicate it contains chemicals that could cause cancer, birth defects or other reproductive harm. Learn more about prop65 at [www.p65warnings.ca.gov](http://www.p65warnings.ca.gov).



10 GLENSHAW STREET  
ORANGEBURG, NY 10962



1.888 DYNAREX



DYNAREX.COM

FIND US



WE CARE LIKE FAMILY™